

=> d his

(FILE 'HOME' ENTERED AT 14:24:42 ON 06 AUG 2001)

FILE 'HCAPLUS' ENTERED AT 14:24:48 ON 06 AUG 2001

L1 73 S ICT
L2 45372 S IMMUNOASSAY? OR IMMUNOCHEMICAL ANALYSI?
L3 191 S IMMUNOCHROMATOG? (L) (ASSAY? OR IMMUNOASSAY?)
L4 114196 S L2 OR ASSAY?
L5 264 S L1 OR L3
L6 101952 S (H OR HAEMOPHILUS)
L7 4047 S (H OR HAEMOPHILUS) (2W) INFLUEN?
L8 1 S L7 AND L5
L9 151595 S ANTIBOD?
L10 181950 S ANTIGEN?
L11 33182 S L9 (L) L10
L12 94 S L11 AND L7
L13 14 S L12 AND L4
L14 6647 S L4 (L) APP?
L15 4 S L13 AND L14
L16 349 S (IMMUNOCHROMAT?)/AB
L17 0 S L16 AND L7
L18 399 S L3 OR L16
L19 77 S L18 AND L9 AND L10
L20 99775 S CARBOHYDRAT?
L21 43853 S POLYSACCHARID?
L22 139446 S L20 OR L21
L23 4400 S L22 (L) L10
L24 5 S L23 AND L19
L25 29421 S STRIP#
L26 0 S L12 AND L25
L27 0 S L7 AND L4 AND L25
L28 2 S L7 AND L4 AND STRIP?/AB
L29 7 SS L8 OR L15 OR L28

FILE 'REGISTRY' ENTERED AT 14:34:33 ON 06 AUG 2001
E GOLD/CN

L30 1 S E3

FILE 'HCAPLUS' ENTERED AT 14:34:40 ON 06 AUG 2001

L31 142553 S L30 OR GOLD OR AU
L32 1159 S L31 AND L4
L33 0 S L32 AND L7
L34 26 S L31 AND L3

FILE 'STNGUIDE' ENTERED AT 14:35:57 ON 06 AUG 2001

FILE 'HCAPLUS' ENTERED AT 14:42:01 ON 06 AUG 2001

L35 0 S L34 AND L7
L36 2 S L34 AND L23
L37 9 S L29 OR L36

=> d .ca 137 1-9

L37 ANSWER 1 OF 9 HCAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 2000:209941 HCAPLUS
 DOCUMENT NUMBER: 132:233987
 TITLE: Process and materials for the rapid detection of
 Streptococcus pneumoniae employing purified
 antigen-specific antibodies
 INVENTOR(S): Moore, Norman James; Fent, Mary Kathleen; Koulchin,
 Vladimir Andrei; Molokova, Elena Valentin
 PATENT ASSIGNEE(S): Binax, Inc., USA
 SOURCE: PCT Int. Appl., 38 pp.
 CODEN: PIXXD2
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 1
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2000016803	A1	20000330	WO 1999-US21505	19990920
W: AT, AU, CA, CH, CN, CZ, DE, DK, ES, FI, GB, HU, IL, IN, JP, KR, LU, MX, NO, NZ, PL, PT, RU, SE, SK, UA, ZA				
RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
AU 9961513	A1	20000410	AU 1999-61513	19990920
EP 1113817	A1	20010711	EP 1999-948305	19990920
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, FI				

PRIORITY APPLN. INFO.: US 1998-156486 A 19980918
 US 1999-397110 A 19990916
 WO 1999-US21505 W 19990920

AB A process is disclosed for obtaining a C-polysaccharide cell wall antigen
 contg. not more than about 10% protein from Streptococcus pneumoniae
 bacteria. The antigen thus obtained is conjugated to a spacer mol., and
 the free end of the latter is then conjugated to a chromatog. affinity
 column. The column is then utilized to purify raw antibodies to S.
 pneumoniae bacteria, thereby producing antigen-specific antibodies. A
 portion of such antibodies is conjugated to a labeling agent which
 displays a visible color change upon reaction of the antibodies with
 their
 antigenic binding partner and embedded in a first zone of an
 immunochromatog. assay device. Another portion of such antibodies is
 bound to the reaction zone of the device which has a view window. When a
 liq. sample, such as patient urine, cerebrospinal fluid or blood is
 applied to the first zone, the conjugate of antibodies and labeling agent
 and the sample move along a flow strip of bibulous material to the
 reaction zone wherein, if the sample contains S. pneumoniae or its cell
 wall antigen, a sandwich is formed among the labeled conjugate, the
 antigen and the bound antibodies and a color change is obsd. The
 immunochromatog. assay thus performed is completed within about 15 min.
 This assay affords a basis for rapid and reliable diagnosis of various
 pathogenic states caused by S. pneumoniae including pneumonia,
 bronchitis,
 otitis media, sinusitis, meningitis, and secondary disease states that
 commonly occur when primary pneumonic infection caused by this bacterium
 persists undiminished over a time period of 3-5 days.

IC ICM A61K039-385
 ICS A61K039-085
 CC 9-10 (Biochemical Methods)

IT **Immunoassay**
 (app., **Immunochromatog. assay** device; process and materials for rapid detection of *Streptococcus pneumoniae* employing purified antigen-specific antibodies)

IT **Polysaccharides**, analysis
 Proteins, general, analysis
 RL: ANT (Analyte); THU (Therapeutic use); ANST (Analytical study); BIOL (Biological study); USES (Uses)
 (process and materials for rapid detection of *Streptococcus pneumoniae* employing purified **antigen**-specific antibodies)

IT **Carbohydrates**, preparation
 RL: ARG (Analytical reagent use); PUR (Purification or recovery); ANST (Analytical study); PREP (Preparation); USES (Uses)
 (process and materials for rapid detection of *Streptococcus pneumoniae* employing purified **antigen**-specific antibodies)

IT 302-01-2, Hydrazine, uses **7440-57-5, Gold**, uses
 RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
 (process and materials for rapid detection of *Streptococcus pneumoniae* employing purified antigen-specific antibodies)

REFERENCE COUNT: 11
 REFERENCE(S): (1) Bennett; The Journal of Immunology 1979, V122(6), P2356 HCAPLUS
 (2) Fischer; European Journal of Biochemistry 1993, P851 HCAPLUS
 (4) Jennings; Biochemistry 1980, V19(20), P4712 HCAPLUS
 (5) Lees; US 5849301 A 1998 HCAPLUS
 (6) Marburg; US 4830852 A 1989 HCAPLUS
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L37 ANSWER 2 OF 9 HCAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 1996:494700 HCAPLUS
 DOCUMENT NUMBER: 125:162734
 TITLE: Methods for detection of an analyte
 INVENTOR(S): Bogart, Gregory R.; Moddel, Garret R.; Maul, Diana M.;
 Etter, Jeffrey B.; Crosby, Mark
 PATENT ASSIGNEE(S): Biostar, Inc., USA
 SOURCE: U.S., 71 pp. Cont.-in-part of U.S. Ser. No. 924343, abandoned.
 CODEN: USXXAM
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 14
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 5541057	A	19960730	US 1993-75952	19930610
AU 9179004	A1	19921021	AU 1991-79004	19910320
AU 653940	B2	19941020		
EP 539383	A1	19930505	EP 1991-910056	19910320
EP 539383	B1	19960918		
R: BE, CH, DE, ES, FR, GB, IT, LI, LU, NL, SE				
JP 05506936	T2	19931007	JP 1991-509344	19910320
ES 2094224	T3	19970116	ES 1991-910056	19910320
US 5639671	A	19970617	US 1995-412600	19950328

US 5629214	A	19970513	US 1995-456040	19950531
US 5869272	A	19990209	US 1995-455652	19950531
JP 10288616	A2	19981027	JP 1998-5911	19980114
JP 2951300	B2	19990920		

PRIORITY APPLN. INFO.:

US 1989-408291	B2	19890918
US 1992-873097	B2	19920424
US 1992-924343	B2	19920731
JP 1990-513789	A3	19900918
EP 1991-910056	A	19910320
WO 1991-US1781	A	19910320
US 1992-923048	B2	19920731
US 1993-75952	A3	19930610
US 1993-76319	B1	19930610

AB This invention relates to devices that produce a detectable attenuation of the spectral characteristic of light impinging on the devices by thin-film phenomena. Interference phenomena are central to the devices and methods of the invention. The presence or amt. of an analyte of interest (e.g., rheumatoid factor, viral antigens, Streptococcus Group A antigen, allergens, HIV I or II, etc.) in a sample (e.g., blood, urine, spinal fluid, gastric wash, vaginal secretions, etc.) is found by using a substrate having an optically active surface exhibiting a first color in response to light impinging thereon and exhibiting a second color comprising a combination of wavelengths of light different from the first color or comprising an intensity of at least one wavelength of light different from the first color in response to the light when the analyte is present on the surface. Then the optically active surface is contacted with a sample potentially comprising the analyte of interest under conditions in which the analyte can interact with the optically active surface to cause the optically active surface to exhibit the second color when the analyte is present. The devices permit detection of extremely small quantities of analyte in a sample, in amts. as low as 0.1 nM, 0.1 ng/mL, 50 fg, or 2 .times. 10³ organisms in a rapid assay that lasts only a few minutes.

IC ICM C12Q001-70
ICS G01N033-53; G01N033-543; G01N021-00

NCL 435005000

CC 9-1 (Biochemical Methods)
Section cross-reference(s): 3, 7, 15, 73

ST optically active surface app biochem analysis; interference film optical app biochem analysis; thin film analyzer body fluid; bacteria detection body fluid app; virus detection body fluid app; **antigen** detection body fluid app; **antibody** detection body fluid app

IT Bacteria
Blood analysis
Body fluid
Cerebrospinal fluid
Chlamydia
Ellipsometers
Escherichia coli
Feces
Films
Immunoassay
Infrared radiation
Interference

Latex
Light
Neisseria meningitidis
Optical detectors
Pericardium
Peritoneum
Pharynx
Pleura
Reflectometers
Respiratory tract
Saliva
Sputum
Stomach
Streptococcus pneumoniae
Ultraviolet radiation
Urine analysis

(app. and methods for anal. using thin-film phenomena)

IT **Haemophilus influenzae**

(type b, app. and methods for anal. using thin-film phenomena)

L37 ANSWER 3 OF 9 HCAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 1996:191556 HCAPLUS

DOCUMENT NUMBER: 124:337288

TITLE: Devices and methods for detection of an analyte based upon light interference

INVENTOR(S): Sandstrom, Torbjorn; Stibler, Lars; Maul; Diana M.

PATENT ASSIGNEE(S): Biostar, Inc., USA

SOURCE: U.S., 69 pp., Cont.-in-part of U. S. Ser. No.

923,268,

abandoned.

CODEN: USXXAM

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 5494829	A	19960227	US 1993-75128	19930610
US 5631171	A	19970520	US 1995-455493	19950531
PRIORITY APPLN. INFO.:			US 1992-923268	19920731
			US 1993-75128	19930610

AB Instrument configured and arranged to detect the presence or amt. of an analyte of interest on the substrate of an optical device is disclosed. The instrument has a source of linearly polarized, monochromatic light positioned at an angle other than Brewster's angle relative to the substrate; and an analyzer positioned at the same angle relative to the substrate at a location suitable for detecting reflected polarized light from the substrate; wherein the analyzer is configured and arranged to approx. maximize the change in intensity of the light reflected from the substrate that is transmitted through the analyzer when a change in mass occurs at the substrate relative to an unreacted surface.

IC ICM G01N033-543

ICS G01N

NCL 436518000

CC 9-1 (Biochemical Methods)

ST app light interference immunoassay; bacteria

antibody antigen detection app
 IT Bacteria
Haemophilus influenzae
 Neisseria meningitidis
 Streptococcus pneumoniae
 (devices and methods for detection of an analyte based upon light interference)
 IT **Immunoassay**
 (enzyme-linked immunosorbent assay, devices and methods for detection of an analyte based upon light interference)

L37 ANSWER 4 OF 9 HCAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 1996:184244 HCAPLUS

DOCUMENT NUMBER: 124:222848

TITLE: Simultaneous detection, identification and differentiation of eubacterial taxa using a hybridization assay

INVENTOR(S): Jannes, Geert; Rossau, Rudi; Van Heuverswyn, Hugo

PATENT ASSIGNEE(S): Innogenetics, N.V., Belg.

SOURCE: PCT Int. Appl., 247 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9600298	A1	19960104	WO 1995-EP2452	19950623
W:	AM, AT, AU, BB, BG, BR, BY, CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GE, HU, IS, JP, KE, KG, KP, KR, KZ, LK, LR, LT, LU, LV, MD, MG, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, TJ, TM, TT			
RW:	KE, MW, SD, SZ, UG, AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG			
CA 2193101	AA	19960104	CA 1995-2193101	19950623
AU 9529246	A1	19960119	AU 1995-29246	19950623
EP 769068	A1	19970423	EP 1995-924923	19950623
R:	AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LI, LU, MC, NL, PT, SE			
BR 9508101	A	19971230	BR 1995-8101	19950623
JP 10501976	T2	19980224	JP 1995-502804	19950623
EP 1088899	A2	20010404	EP 2001-200037	19950623
EP 1088899	A3	20010502		
R:	AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV			
EP 1098006	A1	20010509	EP 2001-200045	19950623
R:	AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV			
EP 1098007	A1	20010509	EP 2001-200046	19950623
R:	AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE			
EP 1091004	A2	20010411	EP 2001-200042	19950624
EP 1091004	A3	20010418		
R:	AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE			

US 6025132	A	20000215	US 1996-765332	19961223
AU 9935035	A1	19990902	AU 1999-35035	19990615
AU 723742	B2	20000907		

PRIORITY APPLN. INFO.:

EP 1994-870106	A	19940624
EP 1995-870032	A	19950407
AU 1995-29246	A3	19950623
EP 1995-924923	A3	19950623
WO 1995-EP2452	W	19950623

AB A method is presented for detection and identification of at least one microorganism, or for the simultaneous detection of several microorganisms

in a sample. The polynucleic acids present in the sample may need to be released, isolated, or concd., and if necessary, the 16S-23S rRNA spacer region, or a part of it, is amplified with at least one suitable primer pair. The polynucleic acids are hybridized with at least one and preferably more than one of the 112 spacer probes or their equiv., under the appropriate hybridization and wash conditions, and/or with a taxon-specific probe derived from any of 103 spacer sequences under the same hybridization and wash conditions. The hybrids formed with each of the probes used under appropriate hybridization and wash conditions are detected, the microorganism(s) present in the sample identified from the differential hybridization signals obtained.

IC ICM C12Q001-68

CC 3-1 (Biochemical Genetics)
Section cross-reference(s): 10

IT Acinetobacter
Acinetobacter baumannii
Bordetella pertussis
Brucella
Campylobacter
Cerebrospinal fluid
Chlamydia
Chlamydia psittaci
Chlamydia trachomatis
Digestive tract
Food analysis
Genitourinary tract
Haemophilus ducreyi
Haemophilus influenzae
Listeria
Listeria monocytogenes
Microorganism
Moraxella catarrhalis
Mycobacterium
Mycobacterium avium
Mycobacterium celatum
Mycobacterium chelonae
Mycobacterium fortuitum
Mycobacterium genavense
Mycobacterium gordonae
Mycobacterium haemophilum
Mycobacterium intracellulare
Mycobacterium kansasii
Mycobacterium malmoense
Mycobacterium marinum
Mycobacterium paratuberculosis
Mycobacterium scrofulaceum

- Mycobacterium simiae
Mycobacterium tuberculosis
Mycobacterium ulcerans
Mycobacterium xenopi
Mycoplasma
Mycoplasma genitalium
Mycoplasma pneumoniae
Neisseria gonorrhoeae
Neisseria meningitidis
Pseudomonas
Pseudomonas aeruginosa
Respiratory tract
Salmonella
Staphylococcus
Staphylococcus aureus
Staphylococcus epidermidis
Streptococcus
Streptococcus agalactiae
Streptococcus pneumoniae
Yersinia enterocolitica
(simultaneous detection, identification and differentiation of
eubacterial taxa using a hybridization assay based on the 16 S-23 S
rRNA spacer region)
- IT 174696-37-8
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(probe CHPS-ICT 1; simultaneous detection, identification and
differentiation of eubacterial taxa using a hybridization assay based
on the 16 S-23 S rRNA spacer region)
- IT 174696-34-5
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(probe CHTR-ICT 1; simultaneous detection, identification and
differentiation of eubacterial taxa using a hybridization assay based
on the 16 S-23 S rRNA spacer region)
- IT 174696-35-6
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(probe CHTR-ICT 2; simultaneous detection, identification and
differentiation of eubacterial taxa using a hybridization assay based
on the 16 S-23 S rRNA spacer region)
- IT 174696-36-7
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(probe CHTR-ICT 3; simultaneous detection, identification and
differentiation of eubacterial taxa using a hybridization assay based
on the 16 S-23 S rRNA spacer region)

L37 ANSWER 5 OF 9 HCAPLUS COPYRIGHT 2001 ACS
ACCESSION NUMBER: 1995:774831 HCAPLUS
DOCUMENT NUMBER: 123:164647
TITLE: Interrupted-flow assay device
INVENTOR(S): Chandler, Howard M.
PATENT ASSIGNEE(S): Smithkline Diagnostics, USA
SOURCE: PCT Int. Appl., 76 pp.
CODEN: PIXXD2
DOCUMENT TYPE: Patent
LANGUAGE: English
FAMILY ACC. NUM. COUNT: 9
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9516208	A1	19950615	WO 1994-US14004	19941206
W: AM, AT, AU, BB, BG, BR, BY, CA, CH, CN, CZ, DE, DK, ES, FI, GB, GE, HU, JP, KE, KG, KP, KR, KZ, LK, LT, LU, LV, MD, MG, MN, MW, NL, NO, NZ, PL, PT, RO, RU, SD, SE, SI, SK, TJ, TT, UA, US, UZ,				
RW: KE, MW, SD, SZ, AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG				
US 5468648	A	19951121	US 1993-163341	19931207
AU 9513008	A1	19950627	AU 1995-13008	19941206
AU 684585	B2	19971218		
EP 733211	A1	19960925	EP 1995-904245	19941206
EP 733211	B1	19980513		
R: BE, CH, DE, ES, FR, GB, IT, LI, NL, SE				
JP 09506177	T2	19970617	JP 1994-516281	19941206
PRIORITY APPLN. INFO.:				
			US 1993-163341	A 19931207
			US 1991-706639	A2 19910529
			US 1992-888831	B2 19920527
			US 1993-40430	A2 19930331
			WO 1994-US14004	W 19941206
AB	The present invention provides chromatog. assay devices that can perform multiple assays simultaneously in the same test strip, as well as methods for their use. One of the assays can be an immunol. assay to detect an antigen, such as human chorionic gonadotropin, while another assay can be a serol. assay to detect an antibody, such as antirubella antibody. An assay device according to the present invention can comprise: (1) a first opposable component including at least one chromatog. medium having a specific binding partner to the first analyte and a specific binding partner to the second analyte immobilized thereto in sep., discrete, non-overlapping zones; and (2) a second opposable component including an absorber. The first and second opposable components are configured such that bringing the first and second opposable components into opposition causes the absorber to come into operable contact with at least one chromatog. medium so that the zone contg. the specific binding partner to the first analyte is functionally divided from the zone contg. the specific binding partner to the second analyte so that both analytes can be detected.			
IC	ICM G01N033-558			
	ICS G01N033-543			
CC	9-1 (Biochemical Methods)			
	Section cross-reference(s): 14, 15			
ST	interrupted flow assay app biochem analysis; antibody antigen detn interrupted flow app; immuno chromatog assay antibody antigen detn; chromatog immuno assay app			
IT	Carbohydrates and Sugars, analysis Glycoproteins, analysis Haptens Hemoglobins Mucoproteins Proteins, analysis RL: ANT (Analyte); ANST (Analytical study) (interrupted-flow assay app. for antibody and antigen detn.)			
IT	7440-57-5, Gold , uses RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses) (interrupted-flow assay app. for antibody and antigen detn.)			

L37 ANSWER 6 OF 9 HCAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 1993:467251 HCAPLUS

DOCUMENT NUMBER: 119:67251

TITLE: Highly sensitive optical immunoassay using enzyme-labeled reagents

INVENTOR(S): Maul, Diana M.; Crider, Debbie G.; Bilodeau, Robert J.; Bogart, Gregory R.

PATENT ASSIGNEE(S): Biostar, Inc., USA

SOURCE: Can. Pat. Appl., 37 pp.

CODEN: CPXXEB

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 14

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
CA 2078897	AA	19930325	CA 1992-2078897	19920923
AU 9179004	A1	19921021	AU 1991-79004	19910320
AU 653940	B2	19941020		
EP 539383	A1	19930505	EP 1991-910056	19910320
EP 539383	B1	19960918		
R: BE, CH, DE, ES, FR, GB, IT, LI, LU, NL, SE				
JP 05506936	T2	19931007	JP 1991-509344	19910320
ES 2094224	T3	19970116	ES 1991-910056	19910320
EP 546222	A1	19930616	EP 1991-308968	19911001
EP 546222	B1	19970910		
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE				
AT 158080	E	19970915	AT 1991-308968	19911001
ES 2109258	T3	19980116	ES 1991-308968	19911001
US 5955377	A	19990921	US 1995-403565	19950417
PRIORITY APPLN. INFO.:			EP 1991-308968	A 19911001
			US 1991-653052	A2 19910211
			EP 1991-910056	A 19910320
			WO 1991-US1781	A 19910320
			US 1992-923090	B2 19920731
			US 1993-75693	B1 19930610

AB A thin-film optical immunoassay device is disclosed which comprises a solid support substrate having an upper and lower surface, an unlabeled ligand antibody layer bound to the substrate, a layer comprising an immobilized enzyme conjugate, complexed with an analyte of interest

and

capable of further interacting with an enzyme-reactive delivery substance to form an insol. reaction product: The enzyme conjugate layer and the unlabeled antibody layer have a measurable increased mass change, the

mass

capable of pptn. by a pptg. agent applied as a substrate thereover. Also provided are a corresponding process for detection of an analyte in a medium and a diagnostic test kit for a thin-film optical immunoassay.

The

invention overcomes the limitations imposed by the prior art use of particulate reagent enhancers. By the use of enzyme-antibody conjugates in place of latex-reagent particles, even more highly sensitive optical thin-film assays can be obtained, particularly with selected substrates for the enzyme which provide insol. pptd. products. The invention

relates

to the use of such enzyme-labeled antibody methods in thin-film assays
for
the detection of low levels of the polysaccharide antigens derived from
the group of bacteria commonly responsible for bacterial infections in
man, e.g. streptococcus. Thus, a conjugate of peroxidase with Igs from
antisera against Neisseria meningitidis (A, C, Y, W135) was dild. in
casein-contg. buffer and mixed with an equal vol. of a diln. of a
cell-free filtrate from a culture of Neisseria meningitidis. The mixt.
was pipeted onto the surface of a silicon wafer already coated with
layers
of silicon nitride, t-polymer siloxane, and purified Ig from the same
anti-N. meningitidis antiserum. TMBlue was used as pptg. substrate; the
ppt. was read by eye and by ellipsometer to confirm the presence of N.
meningitidis. Visually, a 1:20,000 diln. of antigen was clearly resolved
from the neg.; a com. kit's sensitivity cut-off is a 1+ at a 1:4000 diln.
of the same antigen prepn.

IC ICM G01N033-569
CC 9-1 (Biochemical Methods)
ST thin film optical **immunoassay**; enzyme antibody conjugate optical
immunoassay; Neisseria thin film optical **immunoassay**;
Streptococcus thin film optical **immunoassay**
IT Bacteria
Streptococcus pneumoniae
(analyte derived from, detn. of, thin-film optical **immunoassay**
device with enzyme-antibody conjugate for, increased sensitivity in
relation to)

IT **Antigens**
Polysaccharides, analysis
RL: ANST (Analytical study)
(bacteria-derived, detn. of, thin-film optical **immunoassay**
device with enzyme-**antibody** conjugate for, increased
sensitivity in relation to)

IT Ellipsometry
(for thin-film optical **immunoassay** device with
enzyme-antibody conjugate, increased sensitivity in relation to)

IT Immobilization, biochemical
(of unlabeled antibody and enzyme conjugate, for optical thin-film
immunoassay device, increased sensitivity in relation to)

IT Films
(thin, optical **immunoassay** device using, with antibody-enzyme
conjugate, increased sensitivity in relation to)

IT **Immunoassay**
(app., thin-film optical, with antibody-enzyme conjugate,
increased sensitivity in relation to)

IT Immunoglobulins
RL: ANST (Analytical study)
(conjugates, anti-Neisseria meningitidis, with peroxidase, for optical
thin-film **immunoassay** of N. meningitidis antigens)

IT Enzymes
RL: ANST (Analytical study)
(conjugates, with antibodies, for thin-film optical **immunoassay**
device, increased sensitivity in relation to)

IT Antibodies
RL: ANST (Analytical study)
(conjugates, with enzymes, for thin-film optical **immunoassay**
device, increased sensitivity in relation to)

IT Neisseria meningitidis

- Streptococcus
 (group A, analyte derived from, detn. of, thin-film optical
immunoassay device with enzyme-antibody conjugate for,
 increased sensitivity in relation to)
- IT Neisseria meningitidis
 Streptococcus
 (group B, analyte derived from, detn. of, thin-film optical
immunoassay device with enzyme-antibody conjugate for,
 increased sensitivity in relation to)
- IT Neisseria meningitidis
 (group C, analyte derived from, detn. of, thin-film optical
immunoassay device with enzyme-antibody conjugate for,
 increased sensitivity in relation to)
- IT Neisseria meningitidis
 (group W-135, analyte derived from, detn. of, thin-film optical
immunoassay device with enzyme-antibody conjugate for,
 increased sensitivity in relation to)
- IT Neisseria meningitidis
 (group Y, analyte derived from, detn. of, thin-film optical
immunoassay device with enzyme-antibody conjugate for,
 increased sensitivity in relation to)
- IT Antibodies
 RL: ANST (Analytical study)
 (monoclonal, to bacteria, detn. of, optical thin-film
immunoassay device for)
- IT **Haemophilus influenzae**
 (type b, analyte derived from, detn. of, thin-film optical
immunoassay device with enzyme-antibody conjugate for,
 increased sensitivity in relation to)
- IT 54827-17-7, 3,3',5,5'-Tetramethylbenzidine
 RL: ANST (Analytical study)
 (as pptg. agent, in thin-film optical **immunoassay** with
 antibody-enzyme conjugate)
- IT 9003-99-0D, Peroxidase, anti-bacterial antibody conjugates
 RL: ANST (Analytical study)
 (for thin-film optical **immunoassay** device)
- IT 12033-89-5, Silicon nitride, biological studies
 RL: BIOL (Biological study)
 (silicon wafer coated with, optical thin-film **immunoassay** in
 relation to)
- IT 7440-22-4, Silver, biological studies
 RL: BIOL (Biological study)
 (silicon wafer colored with, optical thin-film **immunoassay** in
 relation to)
- IT 7440-21-3, Silicon, biological studies
 RL: BIOL (Biological study)
 (wafer, silicon nitride-coated or silver-colored, optical thin-film
immunoassay in relation to)

L37 ANSWER 7 OF 9 HCAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 1987:64041 HCAPLUS

DOCUMENT NUMBER: 106:64041

TITLE: Solid-phase system and **apparatus** for use in
 ligand-receptor **assays**, particularly

immunoassays

INVENTOR(S): Rubenstein, Albert Samuel

PATENT ASSIGNEE(S): Hybritech, Inc., USA

Hines 09/518,165

SOURCE: Eur. Pat. Appl., 30 pp.
 CODEN: EPXXDW
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 1
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
EP 200381	A1	19861105	EP 1986-302521	19860404
EP 200381	B1	19930127		
R: AT, BE, CH, DE, FR, GB, IT, LI, LU, NL, SE				
CA 1272127	A1	19900731	CA 1986-505788	19860403
FI 8601460	A	19861005	FI 1986-1460	19860404
FI 92257	B	19940630		
FI 92257	C	19941010		
DK 8601538	A	19861005	DK 1986-1538	19860404
DK 171928	B1	19970811		
NO 8601317	A	19861006	NO 1986-1317	19860404
NO 170825	B	19920831		
NO 170825	C	19921209		
AU 8655657	A1	19861016	AU 1986-55657	19860404
AU 605101	B2	19910110		
JP 61292059	A2	19861222	JP 1986-78064	19860404
JP 07113636	B4	19951206		
ES 553724	A1	19890201	ES 1986-553724	19860404
ES 553724	A5	19890331		
EP 437287	A2	19910717	EP 1991-103886	19860404
EP 437287	A3	19910724		
EP 437287	B1	19960717		
R: AT, BE, CH, DE, FR, GB, IT, LI, LU, NL, SE				
AT 85129	E	19930215	AT 1986-302521	19860404
AT 140539	E	19960815	AT 1991-103886	19860404
US 5879881	A	19990309	US 1993-162791	19931207

PRIORITY APPLN. INFO.:

US 1985-720036	19850404
US 1986-847799	19860403
EP 1986-302521	19860404
US 1987-3496	19870115
US 1990-579087	19900907
US 1992-879693	19920505

AB A solid-phase system for use in tittle assays (e.g., immunoassays, nucleic acid probe assays) comprises a porous matrix in which microspheres, bound to a receptor capable of capturing a target ligand, are entrapped. Microspheres with different bound receptors may be entrapped at different sites in the matrix enabling the detn. of .gtoreq.2 analytes in a sample at the same time. Also, internal controls may be incorporated by having zones of microspheres without receptors. For the detection of human chorionic gonadotropin (HCG) in urine, polystyrene microspheres were activated with monoclonal antibodies to HCG and then entrapped in Whatman GF/F glass fiber filters incorporated in an immunoconcn. app. Urine, followed by alk. phosphatase-labeled monoclonal antibodies to HCG, were applied to the filter and drawn through. The filter was washed and indoxyl phosphate was added. A dark blue color developed in the discrete zone of activated microspheres, indicating the presence of HCG.

IC ICM G01N033-549
 ICS G01N033-545; G01N033-551; G01N033-577; G01N033-564; G01N033-569;
 G01N033-74; G01N033-70

CC 9-10 (Biochemical Methods)
 Section cross-reference(s): 2, 15

ST solid phase ligand receptor **assay**; chorionic gonadotropin solid
 phase EIA; microsphere ligand receptor **assay**

IT Bacteria
 Candida albicans
 Chlamydia trachomatis
 Fungi
Haemophilus influenzae
 Neisseria gonorrhoeae
 Parasite
 Trichomonas vaginalis
 Virus, animal
 (detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Receptors
 RL: ANST (Analytical study)
 (matrix-entrapped microsphere-bound, for ligands, detn. of, by
 solid-phase **assay**)

IT Allergens
 Antigens
 RL: ANST (Analytical study)
 (matrix-entrapped microsphere-bound, in solid-phase **immunoassay**
)

IT Latex
 (matrix-entrapped microspheres, receptor-bound, in solid-phase
immunoassay)

IT Blood analysis
 Urine analysis
 (matrix-entrapped receptor-bound microspheres in solid-phase
immunoassay for)

IT Ceramic materials and wares
 Filtering materials
 Glass fibers, uses and miscellaneous
 Polyamide fibers, uses and miscellaneous
 RL: USES (Uses)
 (receptor-bound microspheres entrapped in, for solid-phase
immunoassay)

IT Ligands
 RL: ANST (Analytical study)
 (receptor-bound microspheres for, porous matrix-entrapped, in
 solid-phase **assay**)

IT Antibodies
 RL: PROC (Process)
 (to rubella virus, detn. of, with matrix-entrapped receptor-bound
 microspheres in solid-phase **immunoassay**)

IT Immunoglobulins
 RL: ANT (Analyte); ANST (Analytical study)
 (E, detn. of, matrix-entrapped allergen-bound microspheres in
immunoassay for)

IT Virus, animal
 (adeno-, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Antigens
 RL: ANT (Analyte); ANST (Analytical study)
 (carcinoembryonic, detn. of, matrix-entrapped receptor-bound
 microspheres in **immunoassay** for)

IT Virus, animal
(cytomegalo-, detn. of, matrix-entrapped receptor-bound microspheres
in **immunoassay** for)

IT **Immunochemical analysis**
(enzyme **immunoassay**, solid-phase, matrix-entrapped
antigen- or antibody-bound microspheres in)

IT Streptococcus
(group A, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Streptococcus
(group B, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Virus, animal
(hepatitis, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Virus, animal
(hepatitis A, detn. of, matrix-entrapped receptor-bound microspheres
in **immunoassay** for)

IT Virus, animal
(hepatitis B, detn. of, matrix-entrapped receptor-bound microspheres
in **immunoassay** for)

IT Virus, animal
(hepatitis, non-A, non-B, detn. of, matrix-entrapped receptor-bound
microspheres in **immunoassay** for)

IT Virus, animal
(herpes, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Virus, animal
(human T-cell leukemia, detn. of, matrix-entrapped receptor-bound
microspheres in **immunoassay** for)

IT **Immunochemical analysis**
(**immunoassay**, solid-phase, matrix-entrapped **antigen**
- or **antibody-bound microspheres in**)

IT Virus, animal
(influenza, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Antigens
RL: ANT (Analyte); ANST (Analytical study)
(prostate-specific, detn. of, matrix-entrapped receptor-bound
microspheres in **immunoassay** for)

IT Virus, animal
(respiratory syncytial, detn. of, matrix-entrapped receptor-bound
microspheres in **immunoassay** for)

IT Virus, animal
(rota-, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Virus, animal
(rubella, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT Fetoproteins
RL: ANT (Analyte); ANST (Analytical study)
(.alpha.-, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT 9001-77-8 9002-61-3, Choriogonadotropin 9002-67-9

RL: ANT (Analyte); ANST (Analytical study)
 (detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT 9001-15-4
 RL: ANST (Analytical study)
 (isoenzymes, detn. of, matrix-entrapped receptor-bound microspheres in
immunoassay for)

IT 9002-88-4, Polyethylene 9003-07-0, Polypropylene 9003-53-6,
 Polystyrene
 RL: ANST (Analytical study)
 (matrix-entrapped microspheres, receptor-bound, in solid-phase
immunoassay)

L37 ANSWER 8 OF 9 HCAPLUS COPYRIGHT 2001 ACS
 ACCESSION NUMBER: 1985:75210 HCAPLUS
 DOCUMENT NUMBER: 102:75210
 TITLE: Immunoblot method for identifying surface components,
 determining their cross-reactivity, and investigating
 cell topology: results with **Haemophilus**
influenzae type b
 AUTHOR(S): Loeb, Marilyn R.
 CORPORATE SOURCE: Med. Cent., Univ. Rochester, Rochester, NY, 14642,
 USA
 SOURCE: Anal. Biochem. (1984), 143(1), 196-204
 CODEN: ANBCA2; ISSN: 0003-2697
 DOCUMENT TYPE: Journal
 LANGUAGE: English

AB An outer membrane prep. derived from encapsulated (type b) H. influenzae
 was resolved by SDS-polyacrylamide gel electrophoresis, and the sepd.
 components were then transferred electrophoretically to nitrocellulose.
 The nitrocellulose was cut into vertical **strips**, which were then
 each incubated with rabbit antiserum to the whole bacterium or with the
 same antiserum after absorption with any of the following: the same
 strain
 of H. influenzae, a capsule-deficient mutant of that strain, other
 strains
 of H. influenzae, or other bacteria. The **strips** were then
 incubated with ¹²⁵I-protein A, and the bound antibodies were detected by
 autoradiog. The autoradiograph of the **strip** exposed to
 unadsorbed antisera revealed the identity of those individual outer
 membrane components that bound antibodies. A comparison of the intensity
 of the various bands on this **strip** with those on the
strips exposed to adsorbed antisera was then used to identify: (1)
 surface-exposed components, (2) those components occluded by capsule, and
 (3) cross-reactivity of exposed components. This method should be
 applicable to other cells and subcellular particles. Its major
 disadvantage is that it can provide false neg. results.

CC 9-10 (Biochemical Methods)
 Section cross-reference(s): 10

IT Lipopolysaccharides
 Proteins
 RL: ANT (Analyte); ANST (Analytical study)
 (detection of, of **Haemophilus influenzae** surface
 membrane by gel electrophoresis and immunoblot)

IT **Immunochemical analysis**
 (gel electrophoresis combined with, of lipopolysaccharides and
 proteins)

of **Haemophilus influenzae** surface membrane)
 IT **Haemophilus influenzae**
 (lipopolysaccharides and proteins in surface membranes of, by gel
 electrophoresis and immunoblot)
 IT Electrophoresis and Ionophoresis
 (gel, immunoblot combined with, of lipopolysaccharides and proteins of
Haemophilus influenzae surface membrane)
 IT Cell wall
 (outer membrane, lipopolysaccharides in proteins detection in, of
Haemophilus influenzae surface membrane by gel
 electrophoresis and immunoblotting)

L37 ANSWER 9 OF 9 HCAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 1984:625331 HCAPLUS

DOCUMENT NUMBER: 101:225331

TITLE: Pyridinium 2-azo-p-dimethylaniline chromophore, a
 chromogenic reagent for .beta.-lactamase testing
 compared to nitrocefin

AUTHOR(S): Barlam, T.; Neu, H. C.

CORPORATE SOURCE: Coll. Phys. Surg., Columbia Univ., New York, NY,
 10032, USA

SOURCE: Eur. J. Clin. Microbiol. (1984), 3(3), 185-9

CODEN: EJCMDM; ISSN: 0722-2211

DOCUMENT TYPE: Journal

LANGUAGE: English

AB Pyridinium-2-azo-p-dimethylaniline chromophore was evaluated as a test
 tube, filter paper, and spectrometric assay for detection of
 .beta.-lactamases from gram-pos. and gram-neg. organisms. Although
 useful

for detection of TEM .beta.-lactamases in *Haemophilus influenzae* and
Neisseria gonorrhoeae, it was a poor agent for detecting TEM, OXA, and

PSE enzymes in Enterobacteriaceae. It also proved poor for detecting
 cephalosporinases in *Pseudomonas aeruginosa* and Enterobacteriaceae, and
 penicillinases in *Staphylococcus aureus* when compared to nitrocefin. As

a spectrometric substrate it was equiv. to nitrocefin and cephaloridine
 with

various .beta.-lactamases, but it was not useful in the above assays with
 test tubes or filter paper **strips**.

CC 7-1 (Enzymes)

Section cross-reference(s): 9, 10

IT Bacteria

Haemophilus influenzae

Neisseria gonorrhoeae

(.beta.-lactamase of, detn. of, azodimethylaniline chromophore as
 substrate for)

IT 9001-74-5 9012-26-4

RL: ANT (Analyte); ANST (Analytical study)

(detn. of, azodimethylaniline chromophore as substrate for, of
 bacteria, spectrometric **assay** in relation to)

=> fil wpids

FILE 'WPIDS' ENTERED AT 14:50:46 ON 06 AUG 2001
COPYRIGHT (C) 2001 DERWENT INFORMATION LTD

FILE LAST UPDATED: 01 AUG 2001 <20010801/UP>
MOST RECENT DERWENT UPDATE 200143 <200143/DW>
DERWENT WORLD PATENTS INDEX SUBSCRIBER FILE, COVERS 1963 TO DATE

>>> SDI'S MAY BE RUN ON EVERY UPDATE OR MONTHLY AS OF JUNE 2001.
(EVERY UPDATE IS THE DEFAULT). FOR PRICING INFORMATION
SEE HELP COST <<<

>>> D COST AND SET NOTICE DO NOT REFLECT SUBSCRIBER DISCOUNTS -
SEE HELP COST <<<

>>> FOR UP-TO-DATE INFORMATION ABOUT THE DERWENT CHEMISTRY
RESOURCE, PLEASE VISIT
<http://www.derwent.com/chemistryresource/index.html> <<<

>>> FOR DETAILS OF THE PATENTS COVERED IN CURRENT UPDATES,
SEE <http://www.derwent.com/covcodes.html> <<<

=> d his

(FILE 'WPIDS' ENTERED AT 14:44:55 ON 06 AUG 2001)
DEL HIS Y
L1 110 S IMMUNOCHROMATOG?
L2 33 S ICT
L3 140 S L1 OR L2
L4 24802 S ANTIGEN?
L5 39493 S ANTIBOD?
L6 42 S L3 AND L4 AND L5
L7 1425 S H INFLUENZ? OR HAEMOPHI?
L8 0 S L6 AND L7
L9 0 S L7 AND L3
L10 138 S L4 AND L5 AND L7
L11 27593 S ASSAY? OR IMMUNOASSAY? OR IMMUNOCHEM?
L12 38 S L10 AND L11
L13 1562337 S APP## OR APPARAT? OR STRIP#
L14 1 S L12 AND L13
L15 74565 S AU OR GOLD
L16 2 S L12 AND L15
L17 535 S (CARBOHYDRAT## OR ?SACCHARIDE?) (4A) L4
L18 6 S L12 AND L17
L19 7 S L14 OR L16 OR L18

FILE 'WPIDS' ENTERED AT 14:50:46 ON 06 AUG 2001

=> d .wp 1-7

L19 ANSWER 1 OF 7 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD
AN 2001-328612 [34] WPIDS
DNC C2001-100777

TI Determining the dose response of a human to a bacterial polysaccharide conjugate vaccine comprises measuring an immune response.

DC B04 D16

IN LAFERRIERE, C A J; POOLMAN, J; SLAOUI, M M

PA (SMIK) SMITHKLINE BEECHAM BIOLOGICALS

CYC 94

PI WO 2001030390 A2 20010503 (200134)* EN 35p

RW: AT BE CH CY DE DK EA ES FI FR GB GH GM GR IE IT KE LS LU MC MW MZ
NL OA PT SD SE SL SZ TZ UG ZW

W: AE AG AL AM AT AU AZ BA BB BG BR BY BZ CA CH CN CR CU CZ DE DK DM
DZ EE ES FI GB GD GE GH GM HR HU ID IL IN IS JP KE KG KP KR KZ LC
LK LR LS LT LU LV MA MD MG MK MN MW MX MZ NO NZ PL PT RO RU SD SE
SG SI SK SL TJ TM TR TT TZ UA UG US UZ VN YU ZA ZW

ADT WO 2001030390 A2 WO 2000-EP10733 20001027

PRAI GB 1999-25559 19991028

AB WO 200130390 A UPAB: 20010620

NOVELTY - A method (I) of determining the dose response of a human to a polysaccharide conjugate vaccine, is new,

DETAILED DESCRIPTION - A method (I) of determining the dose response of a human to a polysaccharide conjugate vaccine comprising an immunogenic

carrier protein and a bacterial polysaccharide. The method comprises:

(a) administering a dose of the conjugate vaccine to an infant animal; and

(b) determining the immune response of the animal to the bacterial polysaccharides as a measure of the immune response of a human.

INDEPENDENT CLAIMS are included for the following:

(1) The use of animals in (I); and

(2) A combination vaccine (II) comprising 2 or more pneumococcal capsular **polysaccharide** conjugate **antigens** at an optimal concentration for inducing an optimal anti-polysaccharide **antibody** response when administered to a human;

ACTIVITY - Antibacterial.

MECHANISM OF ACTION - Vaccine.

No supporting data is given.

USE - (II) can be used to treat pneumococcal disease in a human host and for the manufacture of medicaments for treatment (both claimed).

The method (I) is used for testing a vaccine response in an animal model to obtain information on the response to humans to the same vaccine **antigen** (claimed). The animals may be used in potency tests for the lot release of batches of vaccines to ensure that a related response in humans would be acceptable, and in pre-clinical studies to evaluate

the

efficacy of new formulations of conjugate without initially having to conduct human trials. The animal model is used to develop a combination vaccine comprising 2 or more pneumococcal capsular **polysaccharide** conjugate **antigens** at an optimal concentration of inducing an optimal anti-polysaccharide **antibody** response when administered to a human.

ADVANTAGE - The animal model is highly predictive of the human response to the **antigen**.

Dwg.0/12

L19 ANSWER 2 OF 7 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD

AN 2000-442428 [38] WPIDS

DNN N2000-330108 DNC C2000-134585

TI New biosensors for detection of analytes, especially yeast, fungi,

bacteria, or virus, comprises an **antibody**-binding protein layer printed in a pattern onto a polymer film and an **antibody** attached to it.

DC A18 A28 A96 B04 D16 J04 M13 S03
IN EVERHART, D S; KAYLOR, R M; MCGRATH, K
PA (KIMB) KIMBERLY-CLARK WORLDWIDE INC
CYC 89
PI WO 2000036416 A1 20000622 (200038)* EN 34p

RW: AT BE CH CY DE DK EA ES FI FR GB GH GM GR IE IT KE LS LU MC MW NL
OA PT SD SE SL SZ TZ UG ZW

W: AE AL AM AT AU AZ BA BB BG BR BY CA CH CN CR CU CZ DE DK DM EE ES
FI GB GD GE GH GM HR HU ID IL IN IS JP KE KG KP KR KZ LC LK LR LS
LT LU LV MA MD MG MK MN MW MX NO NZ PL PT RO RU SD SE SG SI SK SL
TJ TM TR TT TZ UA UG UZ VN YU ZA ZW

AU 2000017431 A 20000703 (200046)
ADT WO 2000036416 A1 WO 1999-US27727 19991122; AU 2000017431 A AU 2000-17431
19991122
FDT AU 2000017431 A Based on WO 200036416
PRAI US 1998-213713 19981217
AB WO 200036416 A UPAB: 20000811

NOVELTY - A biosensor comprising a polymer film, and an **antibody**-binding protein (ABP) layer printed in a pattern onto the polymer film, where the ABP layer has an **antibody** on it that is specific for an analyte, is new.

DETAILED DESCRIPTION - INDEPENDENT CLAIMS are also included for the following:

(1) making a biosensor comprising printing a pattern of ABP layer with a subsequent layer of **antibody** on the polymer film;

(2) detecting an analyte in a medium, comprising contacting the medium with a biosensing device of the novelty, transmitting light through

the polymer film, and detecting the presence of the analyte bound to the **antibody**, by detecting a pattern form by diffraction of the light; and

(3) a biosensor comprising a polymer film, and an ABP layer printed in a pattern onto the polymer film where the ABP layer is capable of acting as a receptor for an analyte.

USE - The biosensors can be used for the detection of analytes such as bacteria, yeast, fungus, virus, rheumatoid factor, IgG, IgM, IgA and IgE **antibodies**, carcinoembryonic **antigen**, streptococcus Group A **antigen**, viral **antigens**, **antigens** associated with autoimmune disease, allergens, tumor **antigens**, streptococcus Group B **antigen**, human immunodeficiency virus (HIV) I or HIV II **antigen**, **antibodies**, viruses, **antigens** specific to Rous sarcoma virus (RSV) **antigen**, enzyme, hormone, **polysaccharide**, protein, lipids, carbohydrate, drug, nucleic acid, Neisseria meningitides groups A, B, C, Y and W sub 135, Streptococcus pneumoniae, Escherichia coli K1, Haemophilus influenza type B, an **antigen** derived from microorganisms, a hapten, a drug of abuse, a therapeutic drug, an environmental agent, **antigens** specific to hepatitis, or especially Candida sp. or Salmonella sp. (claimed). The biosensors can be used to detect contamination in garments, such as diapers, and to detect contamination by microorganisms. They can also be used on contact lenses, eyeglasses, window panes, pharmaceutical vials, solvent containers, water bottles, or adhesive bandages, to detect contamination.

ADVANTAGE - The method allows a modular production format so that

large rolls of patterned protein may be made for use with different analytes. The final product may be made by exposure to the necessary **antibody**. A diffraction image is produced which can be easily seen with the eye or, optionally, with a sensing device. The biosensors allow detection of extremely small quantities of analyte in a medium in a rapid **assay** lasting only a few minutes. In addition, no signaling or associate electronic components are required in the biosensing devices.
Dwg.0/2

L19 ANSWER 3 OF 7 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD
AN 1999-561863 [47] WPIDS
DNC C1999-163794
TI Composition for inducing a type 2 T cell-independent immune response to an **antigen**, particularly for protection against bacteria.
DC B04 C06 D16
IN MOND, J J; SNAPPER, C M
PA (JACK-N) JACKSON FOUND ADVANCEMENT MILITARY MED
CYC 21
PI WO 9947168 A2 19990923 (199947)* EN 41p
RW: AT BE CH CY DE DK ES FI FR GB GR IE IT LU MC NL PT SE
W: AU CA JP
AU 9930054 A 19991011 (200008)
ADT WO 9947168 A2 WO 1999-US5647 19990315; AU 9930054 A AU 1999-30054 19990315
FDT AU 9930054 A Based on WO 9947168
PRAI US 1998-39247 19980316
AB WO 9947168 A UPAB: 19991116
NOVELTY - Composition (A) comprises a type 2 T cell-independent **antigen** (Ag) conjugated to a lipid (or lipid-containing component) (B) that promotes an immune response to Ag.
DETAILED DESCRIPTION - INDEPENDENT CLAIMS are also included for the following:
(1) pharmaceutical composition or vaccine containing (A) dissolved or suspended in a carrier; and
(2) **assay** system for identifying compositions that stimulate isotype switching and/or a B cell memory response.
ACTIVITY - Antibacterial.
MECHANISM OF ACTION - Induction of a specific immune response.
USE - (A) are used to promote a type 2 T cell-independent immune response to Ag (particularly a bacterial polysaccharide), especially in immunocompromised humans, e.g. patients deficient in T cell, neonates, the elderly, or subjects who are immunodefective because of exposure to viruses or other microorganisms, radiation, cytotoxic chemicals, corticosteroids or other immunosuppressants. (A) are useful in human or veterinary medicine.
ADVANTAGE - (A) elicit a primary response that is about 10 times greater than that produced when Ag and B are administered as separate components. They also induce a secondary response, i.e. a vigorous memory response and/or promotion of isotype switching. (A) are effective even in absence of functional cytokine-producing T cells. Mice were rendered T cell defective by administration of an anti-CD4 **antibody**, then one day later they were injected with 5 micro g of (i) Streptococcus pneumoniae type 14 polysaccharide/lipoprotein D complex or (ii) the polysaccharide only. After 14 days, the immunoglobulin G1 titer (ng/ml)

against the polysaccharide was 5834 in (i) but less than 10 in (ii).
Dwg.0/0

L19 ANSWER 4 OF 7 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD
AN 1999-443884 [37] WPIDS
DNN N1999-331070 DNC C1999-130702
TI Biosensors used to detect analyte, e.g. chemical or biological contamination in garments, e.g. diapers.
DC A18 A23 A26 A96 B04 D16 D22 E13 J04 P75 S03
IN EVERHART, D S; JONES, M L; KAYLOR, R M
PA (KIMB) KIMBERLY-CLARK WORLDWIDE INC
CYC 84
PI WO 9931486 A1 19990624 (199937)* EN 39p
RW: AT BE CH CY DE DK EA ES FI FR GB GH GM GR IE IT KE LS LU MC MW NL
OA PT SD SE SZ UG ZW
W: AL AM AT AU AZ BA BB BG BR BY CA CH CN CU CZ DE DK EE ES FI GB GD
GE GH GM HR HU ID IL IN IS JP KE KG KP KR KZ LC LK LR LS LT LU LV
MD MG MK MN MW MX NO NZ PL PT RO RU SD SE SG SI SK SL TJ TM TR TT
UA UG UZ VN YU ZW
AU 9919205 A 19990705 (199948)
US 6060256 A 20000509 (200030)
EP 1040338 A1 20001004 (200050) EN
R: BE DE ES FR GB IT NL SE
CN 1286753 A 20010307 (200140)
ADT WO 9931486 A1 WO 1998-US26759 19981216; AU 9919205 A AU 1999-19205
19981216; US 6060256 A US 1997-991644 19971216; EP 1040338 A1 EP
1998-963991 19981216, WO 1998-US26759 19981216; CN 1286753 A CN
1998-812255 19981216
FDT AU 9919205 A Based on WO 9931486; EP 1040338 A1 Based on WO 9931486
PRAI US 1997-991644 19971216
AB WO 9931486 A UPAB: 19990914
NOVELTY - Biosensors comprising:
(a) polymer film coated with metal; and
(b) patterned receptor layer printed onto (a) on which is a
receptive
material that specifically binds analyte.
DETAILED DESCRIPTION - INDEPENDENT CLAIMS are also included for:
(1) Methods of detecting analyte using the sensor; and
(2) A method of making the biosensor.
USE - Used to detect analyte. Used particularly in the field of
microcontact printing binders on metal films to produce optical
diffraction biosensors. Used as single tests for detecting analyte or as
multiple test devices. Used for detection of chemical or biological
contamination in garments such as diapers, detection of contamination by
microorganisms in pre-packed foods such as fruit juices and other
beverages, and in health diagnostic applications such as diagnostic kits
for detection of **antigens**, microorganisms, and blood
constituents. May be used on contact lenses, eyeglasses, windowpanes,
pharmaceutical vials, solvent containers, water bottles and plasters to
detect contamination. Used in **immunoassays** for either
antigen or **antibody** detection, for use in direct,
indirect or competitive detection systems, for determination of enzymatic
activity, for detection of small organic molecules (drugs of abuse,
therapeutic drugs, environmental agents) and nucleic acids. Samples
with
patterned **antibody** to *Candida albicans* were prepared by
pretreating **gold**/polyester (10 nm thick) by immersion in 5 mg/ml

phosphate-buffered saline solution (pH 7.2) of beta-casein for 10 minutes.

The sample was rinsed thoroughly with distilled water and dried under a strong nitrogen stream. Contact printing was done using a polydimethylsiloxane stamp with an x,y array of 10 μ m diameter circles. The stamp was coated with a thiolated **antibody** to *Candida albicans* by immersing in a 0.5 mg/ml aqueous solution of **antibody** derivative. After 10 minutes, the stamp was removed and thoroughly dried using a strong stream of nitrogen. Contact printing was done on the casein-treated sample, with exposure times of 1 second to 2 minutes being adequate. After printing, the sample was again rinsed with distilled water

and dried. The sensor sample was exposed to germ tube-bearing cells of *C. albicans* by inoculating tape-stripped adult forearm skin with a concentration of 10⁶ yeast cells/ml and placing the sensor on top of the yeast-containing tape. Transfer of the yeast cells to the sensor was accomplished after only a few seconds of contact. Patterned adhesion of the yeast cells to the sensor was confirmed by microscopic analysis and resulted in a diffraction image upon irradiation with a laser.

ADVANTAGE - Are inexpensive and sensitive devices. Produced by easy, efficient and simple method of contact printing a patterned receptor on

an optically transparent, flexible substrate, that is amenable to continuous processing and does not use self-assembling monolayers. Are simpler than prior art, are not restricted to limitations of self-assembling

monolayers and are easier to manufacture. Are low-cost and disposable and can be mass-produced.

DESCRIPTION OF DRAWING(S) - Schematic representation of metal-plated MYLAR (RTM) film with nutrient backing.

MYLAR film (RTM: polyethylene-terephthalate) 15
metal film 20

receptors specific for microorganism 25
nutrient backing 30

Dwg.1/8

L19 ANSWER 5 OF 7 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD
AN 1997-340536 [31] WPIDS
CR 1991-117618 [16]; 1991-117625 [16]; 1992-300174 [36]; 1992-300183 [36];
1992-349359 [42]; 1994-065810 [08]; 1996-010090 [01]; 1996-076885 [08];
1996-361950 [36]; 1997-288174 [26]; 1999-152762 [13]
DNN N1997-282600 DNC C1997-109335
TI Optimisation of visual signal from optical **assay** device - for detection of analyte, e.g. rheumatoid factor, viral **antigens**, **carbohydrate**, drug or nucleic acid.
DC B04 D16 J04 S02 S03
IN BOGART, G R; ETTER, J B
PA (BIOS-N) BIOSTAR INC
CYC 1
PI US 5639671 A 19970617 (199731)* 69p
ADT US 5639671 A CIP of US 1989-408291 19890918, CIP of US 1992-873097 19920424, CIP of US 1992-923048 19920731, Cont of US 1993-76319 19930610, US 1995-412600 19950328
PRAI US 1993-76319 19930610; US 1989-408291 19890918; US 1992-873097 19920424; US 1992-923048 19920731; US 1995-412600 19950328

AB US 5639671 A UPAB: 19990331

Optimizing a visual signal from an optical **assay** device for the detection of an analyte comprises: (a) providing a substrate having an anti-reflective film on it, with a series of thicknesses (preferably optimal) varied incrementally along the length of the substrate, providing an attachment layer of a chosen thickness on the anti-reflective film and providing a receptive layer of a chosen thickness for the analyte on the attachment layer; (b) contacting the analyte with the receptive layer such that a mass change on the receptive layer results, and (c) determining at least 1 thickness of the series of thicknesses of the anti-reflective film that, in combination with the other layers of the device, produces a visual signal comprising a maximised visual contrast in interference colour upon the change in mass relative to a background interference colour, over a range of concentrations of the analyte.

USE - The method can be used to detect analytes such as rheumatoid factor, immunoglobulin E **antibodies** specific for Birch pollen, carcinoembryonic **antigen**, Streptococcus Group A **antigen**, viral **antigens**, **antigens** associated with autoimmune disease, allergens, a tumour or an infectious microorganism, Streptococcus Group B **antigen**, HIV I or HIV II **antigen**, host response (**antibodies**) to the virus, **antigens** specific to RSV or host response (**antibodies**) to the virus, an **antibody**, **antigen**, enzyme, hormone, **polysaccharide**, protein, lipid, carbohydrate, drug or nucleic acid, analyte derived from causative organisms for meningitis, Neisseria meningitidis groups A, B, C, Y and W135, S. pneumoniae, E. coli K1, **Haemophilus** influenzae type B, **antigen** derived from microorganisms, a hapten, a drug of abuse (including drugs which are unlawful to use without a permit or license), a therapeutic drug, an environmental agent and **antigens** specific to hepatitis.
Dwg.0/18

L19 ANSWER 6 OF 7 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD

AN 1994-200269 [24] WPIDS

DNC C1994-091569

TI Nucleic acid encoding D15 outer membrane protein - esp. of **Haemophilus** influenzae, and related proteins, vectors, antisera etc. useful in vaccines, for diagnosis and for passive immunisation..

DC B04 D16

IN CHONG, P; KLEIN, M; LOOSMORE, S; SIA, D Y C; THOMAS, W; YANG, Y;

LOOSMORE,

S M; SIA, D; YANG, Y P

PA (CONN-N) CONNAUGHT LAB LTD

CYC 28

PI WO 9412641 A1 19940609 (199424)* 161p
RW: AT BE CH DE DK ES FR GB GR IE IT LU MC NL PT SE
W: AU BR CA FI JP KR NO NZ RU UA US

AU 9455565 A 19940622 (199436)

EP 668916 A1 19950830 (199539) EN

R: AT BE CH DE DK ES FR GB GR IE IT LI LU MC NL PT SE

JP 08502417 W 19960319 (199644) 156p

AU 683435 B 19971113 (199803)

BR 9307510 A 19990601 (199927)

JP 2907552 B2 19990621 (199930) 181p
 US 6013514 A 20000111 (200010)
 US 6083743 A 20000704 (200036)
 RU 2141528 C1 19991120 (200041)
 KR 216390 B1 19990816 (200104)
 ADT WO 9412641 A1 WO 1993-CA501 19931123; AU 9455565 A AU 1994-55565
 19931123;
 EP 668916 A1 WO 1993-CA501 19931123, EP 1994-900671 19931123; JP 08502417
 W WO 1993-CA501 19931123, JP 1994-512608 19931123; AU 683435 B AU
 1994-55565 19931123; BR 9307510 A BR 1993-7510 19931123, WO 1993-CA501
 19931123; JP 2907552 B2 WO 1993-CA501 19931123, JP 1994-512608 19931123;
 US 6013514 A WO 1993-CA501 19931123, US 1995-433522 19950912; US 6083743
 A
 Cont of WO 1993-CA501 19931123, Cont of US 1995-433522 19950912, US
 1998-135166 19980818; RU 2141528 C1 WO 1993-CA501 19931123, RU
 1995-117238
 19931123; KR 216390 B1 WO 1993-CA501 19931123, KR 1995-702081 19950523
 FDT AU 9455565 A Based on WO 9412641; EP 668916 A1 Based on WO 9412641; JP
 08502417 W Based on WO 9412641; AU 683435 B Previous Publ. AU 9455565,
 Based on WO 9412641; BR 9307510 A Based on WO 9412641; JP 2907552 B2
 Previous Publ. JP 08502417, Based on WO 9412641; US 6013514 A Based on WO
 9412641; RU 2141528 C1 Based on WO 9412641
 PRAI GB 1992-24584 19921123
 AB WO 9412641 A UPAB: 19940803
 New nucleic acid (I) contains at least a portion coding for a D15 outer
 membrane protein (omp) and has a sequence which is (a) any of 5 (all
 about
 3000bp) reproduced in the specification, or complementary sequences or
 (b)
 hybridisable under stringent conditions with such sequences. Also new are
 (1) recombinant plasmids contg. a segment of (I) at least 18 bp long (and
 opt. expression control elements, (12) proteins (II) encoded by these
 plasmids; (3) purified D15 omp (III); (4) synthetic polypeptides with
 sequences corresp. to (II) or (III), or their variants and mutants which
 retain immunogenicity; (5) antisera or **antibodies** specific for
 (II), (III) or immunologous contg. them; (6) chimeric molecules
 consisting
 of (II) or (III) bonded to another polypeptides, protein or
 polysaccharides.
 USE - (I), (II) and the synthetic polypeptides are useful in
 vaccines
 to protect against **Haemophilus**. D15 can also be used as a
 carrier for **polysaccharide antigens** to form conjugate
 vaccines against other bacteria; to induce immunity to abnormal
 polysaccharides or tumour cells and to generate anti-tumour
antibodies, for coupling to toxins etc. (I), (II) synthetic
 peptides and antisera can also be used diagnostically (in hybridisation
 or
immunoassay procedures) and **antibodies** can be used for
 passive immunisation.
 Dwg.0/11

L19 ANSWER 7 OF 7 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD
 AN 1986-292987 [45] WPIDS
 CR 1991-209899 [29]
 DNN N1986-218851 DNC C1986-126914
 TI Solid phase system for use in ligand-receptor **assays** - e.g. for

detection of human choriogonadotropin **antigen** in urine.

DC A89 B04 D16 J04 S03
 IN RUBENSTEIN, A S
 PA (HYBR-N) HYBRITECH INC
 CYC 18
 PI EP 200381 A 19861105 (198645)* EN 30p
 R: AT BE CH DE FR GB IT LI LU NL SE
 AU 8655657 A 19861016 (198648)
 NO 8601317 A 19861027 (198650)
 FI 8601460 A 19861005 (198703)
 JP 61292059 A 19861222 (198705)
 DK 8601538 A 19861005 (198706)
 ES 8900176 A 19890501 (198924)
 CA 1272127 A 19900731 (199036)
 EP 200381 B1 19930127 (199304) EN 14p
 R: AT BE CH DE FR GB IT LI LU NL SE
 DE 3687589 G 19930311 (199311)
 FI 92257 B 19940630 (199428)
 JP 07113636 B2 19951206 (199602) 13p
 DK 171928 B 19970811 (199739)
 ADT EP 200381 A EP 1986-302521 19860404; JP 61292059 A JP 1986-78064
 19860404;
 ES 8900176 A ES 1986-553724 19860404; EP 200381 B1 EP 1986-302521
 19860404; DE 3687589 G DE 1986-3687589 19860404, EP 1986-302521 19860404;
 FI 92257 B FI 1986-1460 19860404; JP 07113636 B2 JP 1986-78064 19860404;
 DK 171928 B DK 1986-1538 19860404
 FDT DE 3687589 G Based on EP 200381; FI 92257 B Previous Publ. FI 8601460; JP
 07113636 B2 Based on JP 61292059; DK 171928 B Previous Publ. DK 8601538
 PRAI US 1985-720036 19850404
 AB EP 200381 A UPAB: 19940907
 Solid phase system for use in a ligand-receptor **assay** for the
 detection of a selected analyte in a fluid sample comprises a porous
 matrix in which microspheres are entrapped, the microspheres being bound
 with a receptor capable of capturing a target ligand.
 The system is for use in a multiple ligand-receptor **assay**
 for simultaneous detection of at least two selected analytes, distinct
 gps. of microspheres that are entrapped within discrete zones in the
 porous matrix, each gp. of microspheres being bound with a receptor
 capable of capturing a different target ligand.
 The following specific cases are claimed: one gp. of microspheres is
 bound with a receptor capable of capturing a target ligand; and one gp.
 of
 microspheres is bound with either (1) the target ligand or other receptor
 substance as a positive control; or (2) a substance incapable of
 capturing
 the target ligand without a bound component, as a negative control for
 the
 detection the analyte. **Appts.** is also provided.
 USE/ADVANTAGE - The system is useful for **assay** of
antigens such as prostatic acid phosphatase, prostate-specific
antigen alpha foetoprotein, carcinoembryonic **antigen**,
 lentinising hormone and creatine kinase isoenzyme. It is also useful for
 the **assay** of viruses, bacteria, parasites or fungi or associated
antigens or **antibodies**, e.g. rubella, rota, hepatitis,
 herpes and influenza viruses, HTLV, group A and group B streptococcus,
 Neisseria gonorrhea, Trichomonas vaginalis, Candida albicans, Chamydia
 trachomatis and **Haemophilus** influenza. The system may also be

Hines 09/518,165

employed in nucleic acid probe technology.

0/6

Dwg.0/6

<-----User Break----->

=> fil medline

FILE 'MEDLINE' ENTERED AT 14:58:00 ON 06 AUG 2001

FILE LAST UPDATED: 30 JUL 2001 (20010730/UP). FILE COVERS 1958 TO DATE.

On April 22, 2001, MEDLINE was reloaded. See HELP RLOAD for details.

MEDLINE now contains new records from the former NLM HEALTH STAR database. These records have an Entry Date and Update Date of 20010223.

MEDLINE thesauri in the /CN, /CT, and /MN fields incorporate the MeSH 2001 vocabulary. Enter HELP THESAURUS for details.

The OLDMEDLINE file segment now contains data from 1958 through 1965. Enter HELP CONTENT for details.

Left, right, and simultaneous left and right truncation are available in the Basic Index. See HELP SFIELDS for details.

THIS FILE CONTAINS CAS REGISTRY NUMBERS FOR EASY AND ACCURATE SUBSTANCE IDENTIFICATION.

=> d his

```
(FILE 'MEDLINE' ENTERED AT 14:51:27 ON 06 AUG 2001)
      DEL HIS Y
L1      226 S IMMUNOCHROMATOGRAP?
      E HAEMOPHILUS/CT
L2      2001 S HAEMOPHILUS/CT
L3      0 S L1 AND L2
      E IMMUNOCHEMISTRY/CT
      E E3+AL
      E E3+ALL
L4      258309 S IMMUNOCHEMISTRY+NT/CT
      E CHROMATOGRAPHY/CT
      E E3+ALLK
      E E3+ALL
L5      0 S CHROMATOGRAPY/CT
L6      26123 S CHROMATOGRAPHY/CT
L7      7 S L6 AND L2
L8      0 S L7 AND L4
L9      10860 S ASSAY/CT OR IMMUNOASSAY/CT
L10     10860 S ASSAYS/CT OR IMMUNOASSAY/CT
L11     0 S L10 AND L7
L12     2 S L7 AND (ASSAY? OR IMMUNOASSAY?)
L13     99901 S STRIP# OR APP## OR APPARATUS
L14     12 S L13 AND L2
L15     0 S L14 AND (L6 OR L10)
      E GOLD/CT
L16     6189 S GOLD/CT
L17     1 S L16 AND L2
L18     15 S L12 OR L14 OR L17
```

FILE 'MEDLINE' ENTERED AT 14:58:00 ON 06 AUG 2001

=> d .med 118 1-15

L18 ANSWER 1 OF 15 MEDLINE
 AN 2000087378 MEDLINE
 DN 20087378 PubMed ID: 10618273
 TI Immunoglobulins in nasal secretions of healthy humans: structural integrity of secretory immunoglobulin A1 (IgA1) and occurrence of neutralizing antibodies to IgA1 proteases of nasal bacteria.
 AU Kirkeby L; Rasmussen T T; Reinholdt J; Kilian M
 CS Department of Medical Microbiology and Immunology, University of Aarhus, DK-8000 Aarhus C, Denmark.
 SO CLINICAL AND DIAGNOSTIC LABORATORY IMMUNOLOGY, (2000 Jan) 7 (1) 31-9. Journal code: CB7; 9421292. ISSN: 1071-412X.
 CY United States
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 200003
 ED Entered STN: 20000314
 Last Updated on STN: 20000314
 Entered Medline: 20000302
 AB Certain bacteria, including overt pathogens as well as commensals, produce immunoglobulin A1 (IgA1) proteases. By cleaving IgA1, including secretory IgA1, in the hinge region, these enzymes may interfere with the barrier functions of mucosal IgA antibodies, as indicated by experiments in vitro. Previous studies have suggested that cleavage of IgA1 in nasal secretions may be associated with the development and perpetuation of atopic disease. To clarify the potential effect of IgA1 protease-producing bacteria in the nasal cavity, we have analyzed immunoglobulin isotypes in nasal secretions of 11 healthy humans, with a focus on IgA, and at the same time have characterized and quantified IgA1 protease-producing bacteria in the nasal flora of the subjects. Samples in the form of nasal wash were collected by using a washing liquid that contained lithium as an internal reference. Dilution factors and, subsequently, concentrations in undiluted secretions could thereby be calculated. IgA, mainly in the secretory form, was found by enzyme-linked immunosorbent **assay** to be the dominant isotype in all subjects, and the vast majority of IgA (median, 91%) was of the A1 subclass, corroborating results of previous analyses at the level of immunoglobulin-producing cells. Levels of serum-type immunoglobulins were low, except for four subjects in whom levels of IgG corresponded to 20 to 66% of total IgA. Cumulative levels of IgA, IgG, and IgM in undiluted secretions ranged from 260 to 2,494 (median, 777) &mgr;g ml⁻¹. IgA1 protease-producing bacteria (*Haemophilus influenzae*, *Streptococcus pneumoniae*, or *Streptococcus mitis* biovar 1) were isolated from the nasal cavities of seven subjects at 2.1×10^3 to 7.2×10^6 CFU per ml of undiluted secretion, corresponding to 0.2 to 99.6% of the flora. Nevertheless, alpha-chain fragments characteristic of IgA1 protease

activity were not detected in secretions from any subject by immunoblotting. Neutralizing antibodies to IgA1 proteases of autologous isolates were detected in secretions from five of the seven subjects but not in those from two subjects harboring IgA1 protease-producing *S. mitis* biovar 1. alpha-chain fragments different from Fc(alpha) and Fd(alpha) were detected in some samples, possibly reflecting nonspecific

proteolytic

activity of microbial or host origin. These results add to previous evidence for a role of secretory immunity in the defense of the nasal mucosa but do not help identify conditions under which bacterial IgA1 proteases may interfere with this defense.

CT Check Tags: Human; Support, Non-U.S. Gov't

Adult

*Antibodies: IM, immunology

Bacteria: EN, enzymology

*Bacteria: IM, immunology

Child

Chromatography

Haemophilus: ME, metabolism

*IgA: CH, chemistry

*IgA, Secretory: CH, chemistry

Immunoblotting

*Immunoglobulins: AN, analysis

Middle Age

*Nasal Mucosa: IM, immunology

Neutralization Tests

Protease Inhibitors: PD, pharmacology

*Serine Endopeptidases: IM, immunology

Streptococcus: ME, metabolism

L18 ANSWER 2 OF 15 MEDLINE

AN 1998395109 MEDLINE

DN 98395109 PubMed ID: 9727006

TI Purification and characterization of wild-type and mutant "classical" nitroreductases of *Salmonella typhimurium*. L33R mutation greatly diminishes binding of FMN to the nitroreductase of *S. typhimurium*.

AU Watanabe M; Nishino T; Takio K; Sofuni T; Nohmi T

CS Division of Genetics and Mutagenesis, National Institute of Health Sciences, 1-18-1 Kamiyoga, Setagaya-ku, Tokyo 158-8501, Japan.

SO JOURNAL OF BIOLOGICAL CHEMISTRY, (1998 Sep 11) 273 (37) 23922-8.

Journal code: HIV; 2985121R. ISSN: 0021-9258.

CY United States

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Priority Journals

EM 199810

ED Entered STN: 19981021

Last Updated on STN: 19981021

Entered Medline: 19981013

AB "Classical" nitroreductase of *Salmonella typhimurium* is a flavoprotein that catalyzes the reduction of nitroaromatics to metabolites that are toxic, mutagenic, or carcinogenic. This enzyme represents a new class of flavin-dependent enzymes, which includes nitroreductases of *Enterobacter cloacae* and *Escherichia coli*, flavin oxidoreductase of *Vibrio fischeri*, and NADH oxidase of *Thermus thermophilus*. To investigate the structure-function relation of this class of enzymes, the gene encoding a mutant nitroreductase was cloned from *S. typhimurium* strain TA1538NR, and

the enzymatic properties were compared with those of the wild-type. DNA sequence analysis revealed a T to G mutation in the mutant nitroreductase gene, predicting a replacement of leucine 33 with arginine. In contrast to the wild-type enzyme, the purified protein with a mutation of leucine 33 to arginine has no detectable nitroreductase activities in the standard assay conditions and easily lost FMN by dialysis or ultrafiltration. In the presence of an excess amount of FMN, however, the mutant protein exhibited a weak but measurable enzyme activity, and the substrate specificity was similar to that of the wild-type enzyme. Possible mechanisms by which the mutation greatly diminishes binding of FMN to the nitroreductase are discussed.

CT Check Tags: Comparative Study; Support, Non-U.S. Gov't

Amino Acid Sequence

Base Sequence

Binding Sites

Chromatography

Chromatography, Gel

Chromatography, Ion Exchange

Cloning, Molecular

Durapatite

Enterobacter cloacae: EN, enzymology

Escherichia coli: EN, enzymology

*FMN: ME, metabolism

Haemophilus: EN, enzymology

Kinetics

Molecular Sequence Data

Molecular Weight

Mycoplasma: EN, enzymology

*Nitroreductases: GE, genetics

Nitroreductases: IP, isolation & purification

*Nitroreductases: ME, metabolism

*Point Mutation

*Salmonella typhimurium: EN, enzymology

Salmonella typhimurium: GE, genetics

Sequence Alignment

Sequence Homology, Amino Acid

Thermus thermophilus: EN, enzymology

Vibrio: EN, enzymology

L18 ANSWER 3 OF 15 MEDLINE

AN 95298926 MEDLINE

DN 95298926 PubMed ID: 7779956

TI Evaluation of a commercial system for the identification of gram-negative,

nonfermenting bacteria of veterinary importance.

AU Salmon S A; Watts J L; Walker R D; Yancey R J Jr

CS Animal Health Discovery Research, Upjohn Company, Kalamazoo, MI 49001, USA.

SO JOURNAL OF VETERINARY DIAGNOSTIC INVESTIGATION, (1995 Jan) 7 (1) 161-4. Journal code: A2D; 9011490. ISSN: 1040-6387.

CY United States

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Priority Journals

EM 199507

ED Entered STN: 19950726

Last Updated on STN: 19950726

Entered Medline: 19950718

CT Check Tags: Animal

*Bacteriological Techniques: VE, veterinary

*Gram-Negative Bacteria: CL, classification

*Gram-Negative Bacteria: IP, isolation & purification

*Gram-Negative Facultatively Anaerobic Rods: CL, classification

*Gram-Negative Facultatively Anaerobic Rods: IP, isolation & purification

Haemophilus

Reagent Strips

L18 ANSWER 4 OF 15 MEDLINE

AN 92112340 MEDLINE

DN 92112340 PubMed ID: 1730505

TI Clustering of an outer membrane adhesin of *Haemophilus parainfluenzae*.

AU Liljemark W F; Bloomquist C G; Lai C H

CS Department of Diagnostic and Surgical Sciences, School of Dentistry,
University of Minnesota, Minneapolis 55455-0329.

NC DE07014 (NIDCR)

R37-DE04614 (NIDCR)

SO INFECTION AND IMMUNITY, (1992 Feb) 60 (2) 687-9.

Journal code: G07; 0246127. ISSN: 0019-9567.

CY United States

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Priority Journals

EM 199202

ED Entered STN: 19920308

Last Updated on STN: 20000303

Entered Medline: 19920218

AB *Haemophilus parainfluenzae* synthesizes an outer membrane protein adhesin which mediates binding to oral streptococci, salivary pellicle, and neuraminidase-treated erythrocytes. An indirect gold labeling technique and immunoelectron microscopy verified the location of this outer membrane

protein. Further, a clustering of gold particles was observed in irregular patches at the cell surface.

CT Check Tags: Support, U.S. Gov't, P.H.S.

*Bacterial Adhesion

*Bacterial Outer Membrane Proteins: AN, analysis

Gold

***Haemophilus: CH, chemistry**

L18 ANSWER 5 OF 15 MEDLINE

AN 91263343 MEDLINE

DN 91263343 PubMed ID: 2048282

TI Improved protection of swine from pleuropneumonia by vaccination with proteinase K-treated outer membrane of *Actinobacillus* (*Haemophilus*) pleuropneumoniae.

AU Chiang Y W; Young T F; Rapp-Gabrielson V J; Ross R F

CS Veterinary Medical Research Institute, College of Veterinary Medicine,
Iowa State University, Ames 50011.

SO VETERINARY MICROBIOLOGY, (1991 Mar) 27 (1) 49-62.

Journal code: XBW; 7705469. ISSN: 0378-1135.

CY Netherlands

DT Journal; Article; (JOURNAL ARTICLE)

LA English
 FS Priority Journals
 EM 199107
 ED Entered STN: 19910802
 Last Updated on STN: 20000303
 Entered Medline: 19910718

AB The immunogenic and protective potentials of an outer membrane-enriched fraction (OM) from a serotype 5 strain of Actinobacillus (Haemophilus) pleuropneumoniae (APP) and the same OM degraded with proteinase K or periodate were evaluated in swine. Groups of pigs were vaccinated with two doses of OM, proteinase K-treated OM (P-OM), periodate-treated OM (PI-OM), or placebo vaccine and challenged intranasally with the homologous strain of APP. Results from triplicate experiments indicated that proteinase K treatment of OM resulted in an improved efficacy. This improved efficacy of P-OM vaccine over untreated OM vaccine was evidenced not only by less severe lung lesions in P-OM vaccinated pigs but also by significant reduction (P less than 0.05) in the number of P-OM vaccinated pigs which developed lung lesions upon challenge with APP. Assessment of sera from vaccinated animals by immunoblotting, complement fixation test, or ELISA indicated that the immunogenicity of some but not all protein or carbohydrate components were reduced (or eliminated) by proteinase K and periodate treatments respectively.

CT Check Tags: Animal; Female; Male; Support, Non-U.S. Gov't
 *Actinobacillus: IM, immunology
 Actinobacillus Infections: PC, prevention & control
 *Actinobacillus Infections: VE, veterinary
 Antibodies, Bacterial: BI, biosynthesis
 Antigens, Bacterial: IM, immunology
 Bacterial Outer Membrane Proteins: IM, immunology
 Complement Fixation Tests
 Electrophoresis, Polyacrylamide Gel
 Endopeptidase K
 Enzyme-Linked Immunosorbent Assay
 Haemophilus: IM, immunology
 Haemophilus Infections: PC, prevention & control
 Haemophilus Infections: VE, veterinary
 Pleuropneumonia: PC, prevention & control
 *Pleuropneumonia: VE, veterinary
 Random Allocation
 Serine Endopeptidases: ME, metabolism
 Swine
 *Swine Diseases: PC, prevention & control
 *Vaccination: VE, veterinary

L18 ANSWER 6 OF 15 MEDLINE
 AN 88047453 MEDLINE
 DN 88047453 PubMed ID: 3674414
 TI Preparation of cell envelopes of large numbers of individual bacterial strains with the use of an automatic cell disruptor.
 AU Van Alphen L; Romijn C; Brandt H; Geelen L; Zanen H C
 CS Department of Medical Microbiology, Academic Medical Centre, Amsterdam, The Netherlands.
 SO ANALYTICAL BIOCHEMISTRY, (1987 Oct) 166 (1) 36-40.

Journal code: 4NK; 0370535. ISSN: 0003-2697.

CY United States
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 198711
 ED Entered STN: 19900305
 Last Updated on STN: 19900305
 Entered Medline: 19871124

AB Analysis of the cell envelopes of large numbers of bacterial strains is used for the epidemiological and taxonomic investigation of clinical, veterinarian, and ecological isolates. Isolation of cell envelopes requires lysis of the bacteria. We developed an **apparatus** to disrupt bacterial cells of 200 different isolates in suspension by ultrasonication automatically. It is composed of modified standard laboratory equipment (fraction collector, cooling unit, pump), a standard ultrasonifier, and a newly designed control unit, which includes a sampler. This **apparatus** was applied to the analysis of cell envelope proteins of 96 Haemophilus influenzae strains on sodium dodecyl sulfate-polyacrylamide gel electrophoresis within 3 days after the first culture.

CT *Bacterial Outer Membrane Proteins: AN, analysis
 Chemistry, Analytical: IS, instrumentation
 Electrophoresis, Polyacrylamide Gel
Haemophilus: AN, analysis

L18 ANSWER 7 OF 15 MEDLINE
 AN 87132255 MEDLINE
 DN 87132255 PubMed ID: 3493113
 TI Rapid determination of X/V growth requirements of Haemophilus species in broth.
 AU Inzana T J; Clarridge J; Williams R P
 NC AI07145 (NIAID)
 SO DIAGNOSTIC MICROBIOLOGY AND INFECTIOUS DISEASE, (1987 Feb) 6 (2) 93-100.

Journal code: DMI; 8305899. ISSN: 0732-8893.

CY United States
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 198704
 ED Entered STN: 19900303
 Last Updated on STN: 19970203
 Entered Medline: 19870423

AB A broth system was developed for rapid identification of the requirement for X factor (hemin), or V factor (NAD), or both for growth of Haemophilus species. This system was compared to growth around paper discs/**strips** impregnated with factors X and/or V. The broth system consisted of three tubes, each containing brain-heart infusion broth supplemented with V factor, X factor, or both. Each tube was inoculated with a saline suspension of an Haemophilus isolate, and the broths were shaken for aeration at 37 degrees C. Under these conditions turbidity or clumping was usually evident after 4-5 hr only in the broth(s) containing the required supplement(s). A few strains requiring only V factor required overnight incubation. One hundred fifty-six Haemophilus isolates were

tested for growth around supplemented discs/**strips** or in supplemented broths: 129 were H. influenzae/aegypticus, 25 were of various species that required only V factor, and 2 were H. aphrophilus. Ten of 89 H. influenzae isolates from the respiratory tract were misidentified by satellitism. All isolates were correctly identified by growth in supplemented broths. The cost of the broth assay was about 60 cents/test, whereas the satellite assay cost about 120 cents/test. Serotyping and antibiotic sensitivity testing could be performed directly from the broth culture. Determination of X and/or V requirement by Haemophilus species with supplemented broths was sensitive, rapid, and inexpensive.

CT Check Tags: Animal; Comparative Study; Human; Support, U.S. Gov't, P.H.S.
 Blood: MI, microbiology
 Cerebrospinal Fluid: MI, microbiology
 Costs and Cost Analysis
 Culture Media
 Ear: MI, microbiology
 Eye: MI, microbiology
 *Haemophilus: GD, growth & development
 Haemophilus: IP, isolation & purification
 *Haemophilus influenzae: GD, growth & development
 Haemophilus influenzae: IP, isolation & purification
 *Heme: AA, analogs & derivatives
 *Hemin: PD, pharmacology
 Kinetics
 *NAD: PD, pharmacology
 Respiratory System: MI, microbiology

L18 ANSWER 8 OF 15 MEDLINE
 AN 86168818 MEDLINE
 DN 86168818 PubMed ID: 3514664
 TI Comparison of a new commercially prepared porphyrin test and the conventional satellite test for the identification of Haemophilus species that require the X factor.
 AU Gadberry J L; Amos M A
 SO JOURNAL OF CLINICAL MICROBIOLOGY, (1986 Mar) 23 (3) 637-9.
 Journal code: HSH; 7505564. ISSN: 0095-1137.
 CY United States
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 198605
 ED Entered STN: 19900321
 Last Updated on STN: 19900321
 Entered Medline: 19860514
 AB A test with a commercially developed porphyrin test agar was examined for the identification of Haemophilus spp. The porphyrin test agar method was compared with the conventional paper **strip** satellite method in tests with 187 isolates and was found to be easier to perform and interpret, giving a sensitivity of 98.7% and specificity of 94.7%.
 CT Check Tags: Comparative Study
 Agar
 *Bacteriological Techniques
 Culture Media
 *Haemophilus: CL, classification
 Haemophilus: GD, growth & development
 Haemophilus: ME, metabolism

*Heme: AA, analogs & derivatives
 *Hemin: PD, pharmacology
 *Porphyrins: BI, biosynthesis

L18 ANSWER 9 OF 15 MEDLINE
 AN 85073854 MEDLINE
 DN 85073854 PubMed ID: 6509390
 TI Porcine haemophili and actinobacilli: characterization by means of API test **strips** and possible taxonomic implications.
 AU O'Reilly T; Rosendal S; Niven D F
 SO CANADIAN JOURNAL OF MICROBIOLOGY, (1984 Oct) 30 (10) 1229-38.
 Journal code: CJ3; 0372707. ISSN: 0008-4166.
 CY Canada
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 198502
 ED Entered STN: 19900320
 Last Updated on STN: 19900320
 Entered Medline: 19850221
 AB Thirty Haemophilus strains and six Actinobacillus strains, all of porcine origin, were examined for their biochemical reactivity on API 20E and API ZYM test **strips** using dense cell suspensions (supplemented with NAD as appropriate) as **strip** inocula. When combined with a test for V-factor dependency, the use of both **strips** allowed adequate differentiation of closely related organisms. Numerical taxonomic analysis of the data demonstrated that the majority of the haemophili and actinobacilli studied could be placed in one of four major clusters; these clusters contained, respectively, the H. pleuropneumoniae--A. pleuropneumoniae strains, the H. parasuis strains, strains belonging to Haemophilus taxon "minor group," and strains belonging to an unusual group of mannitol-positive, urease-negative haemophili. A representative of Haemophilus species taxon C and an unusual Actinobacillus isolate appeared to be comparatively unrelated to organisms in the four major clusters. Although it may, on occasion, be difficult to place an unusual isolate in any one particular group, owing to the uncertain taxonomy of some of these organisms, it is concluded that API test **strips** can serve as useful tools for the characterization and differentiation of porcine haemophili and actinobacilli.
 CT Check Tags: Animal; Support, Non-U.S. Gov't
 *Actinobacillus: CL, classification
 Actinobacillus: EN, enzymology
 *Haemophilus: CL, classification
 Haemophilus: EN, enzymology
 *Indicators and Reagents
 *Reagent Strips
 *Swine: MI, microbiology

L18 ANSWER 10 OF 15 MEDLINE
 AN 84207971 MEDLINE
 DN 84207971 PubMed ID: 6373253
 TI A modified method for differentiation of Haemophilus influenzae from

Haemophilus parainfluenzae.

AU Santanam P
 SO EUROPEAN JOURNAL OF CLINICAL MICROBIOLOGY, (1984 Apr) 3 (2) 150-1.
 Journal code: EMY; 8219582. ISSN: 0722-2211.
 CY GERMANY, WEST: Germany, Federal Republic of
 DT Letter
 LA English
 FS Priority Journals
 EM 198407
 ED Entered STN: 19900320
 Last Updated on STN: 19900320
 Entered Medline: 19840716
 CT *Bacteriological Techniques
 *Haemophilus: IP, isolation & purification
 *Haemophilus influenzae: IP, isolation & purification
 *Indicators and Reagents
 *Reagent Strips

L18 ANSWER 11 OF 15 MEDLINE
 AN 83298942 MEDLINE
 DN 83298942 PubMed ID: 6412201
 TI An easily prepared and accurate test **strip** for the detection of
 beta-lactamase production by both gram negative and positive organisms.
 AU Mugg P A
 SO PATHOLOGY, (1983 Apr) 15 (2) 175-6.
 Journal code: OTA; 0175411. ISSN: 0031-3025.
 CY Australia
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 198310
 ED Entered STN: 19900319
 Last Updated on STN: 19900319
 Entered Medline: 19831021
 AB A simple method is described for the preparation of a paper test
strip for the rapid detection of the beta-lactamase enzymes of
 Haemophilus sp., N. gonorrhoeae, S. aureus and S. epidermidis. The test
strips were compared with the chromogenic cephalosporin technique
 for beta-lactamase detection and found to give 100% agreement. The
 beta-lactamase detection **strips** were easily prepared, very
 inexpensive, required no special equipment and could be stored for up to
 6
 mth.
 CT **Haemophilus: EN, enzymology**
 *Indicators and Reagents
 Neisseria gonorrhoeae: EN, enzymology
 *Reagent Strips
 Staphylococcus: EN, enzymology
 Staphylococcus aureus: EN, enzymology
 *beta-Lactamases: AN, analysis

L18 ANSWER 12 OF 15 MEDLINE
 AN 83109707 MEDLINE
 DN 83109707 PubMed ID: 6337192
 TI Biotyping of Haemophilus using API 10S--an epidemiological tool?.
 AU Mehtar S; Afshar S A
 SO JOURNAL OF CLINICAL PATHOLOGY, (1983 Jan) 36 (1) 96-9.

Journal code: HT3; 0376601. ISSN: 0021-9746.

CY ENGLAND: United Kingdom

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Abridged Index Medicus Journals; Priority Journals

EM 198303

ED Entered STN: 19900318
Last Updated on STN: 19900318
Entered Medline: 19830311

CT Check Tags: Comparative Study; Human
Adolescence
Adult
Age Factors
Bacteriological Techniques
Child
Child, Preschool

***Haemophilus: CL, classification**
Haemophilus: IP, isolation & purification
Haemophilus: ME, metabolism
Haemophilus Infections: MI, microbiology
Haemophilus influenzae: CL, classification
Haemophilus influenzae: IP, isolation & purification
Haemophilus influenzae: ME, metabolism

***Indicators and Reagents**

***Reagent Strips**

L18 ANSWER 13 OF 15 MEDLINE

AN 82167977 MEDLINE

DN 82167977 PubMed ID: 7040443

TI Evaluation of the rapid penicillinase paper **strip** test for detection of beta-lactamase.

AU Oberhofer T R; Towle D W

SO JOURNAL OF CLINICAL MICROBIOLOGY, (1982 Feb) 15 (2) 196-9.
Journal code: HSH; 7505564. ISSN: 0095-1137.

CY United States

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Priority Journals

EM 198206

ED Entered STN: 19900317
Last Updated on STN: 19900317
Entered Medline: 19820621

AB The penicillin-starch paper **strip** method was compared with the acidometric and iodometric methods for assaying beta-lactamase production,
using fresh isolates of clinically important bacteria. Results obtained by
the three methods were compared for rapidity, accuracy, and stability of reagents. Of the 210 isolates tested by the paper **strip** method, 301 isolates tested by the acidometric method, and 117 isolates tested by the iodometric method, all were in perfect agreement with the disk diffusion susceptibility test except one strain each of Haemophilus influenzae, Staphylococcus aureus, and Staphylococcus epidermidis. The H. influenzae isolate was penicillin resistant and failed to give a positive test for beta-lactamase in all three tests. The staphylococci (intermediate and resistant in susceptibility, respectively) failed to give a positive test for beta-lactamase with the iodometric method. The

results of the paper **strip** method, in which 3,241 strains representing nine species of bacteria were used, correlated completely with disk susceptibility tests except for 2 and 69 strains, respectively, of penicillin-resistant, beta-lactamase-negative *H. influenzae* and *H. parainfluenzae*. The results of this study indicate that the paper **strip** method is accurate, simple to perform, extremely economical, and uses materials that are stable when stored frozen. It is eminently suitable for routine laboratory use.

CT Check Tags: Human; Support, Non-U.S. Gov't

*Bacteria: EN, enzymology

*Bacteriological Techniques

Haemophilus: EN, enzymology

Penicillinase: AN, analysis

Reagent Strips

Staphylococcus: EN, enzymology

*beta-Lactamases: AN, analysis

L18 ANSWER 14 OF 15 MEDLINE

AN 77238680 MEDLINE

DN 77238680 PubMed ID: 888594

TI [Polyacrylamid-gel-electrophoresis of *Haemophilus* proteins (author's transl)].

Elektrophoretische Auftrennung von *Haemophilus*-Proteinen im Polyacrylamid-Gel.

AU Neumann U; Hinz K H

SO ZENTRALBLATT FUR BAKTERIOLOGIE, PARASITENKUNDE, INFESTIONSKRANKHEITEN UND HYGIENE. ERSTE ABTEILUNG ORIGINALE. REIHE A: MEDIZINISCHE MIKROBIOLOGIE UND PARASITOLOGIE, (1977 Jun) 238 (2) 244-50.
Journal code: Y52; 0331570. ISSN: 0300-9688.

CY GERMANY, WEST: Germany, Federal Republic of

DT Journal; Article; (JOURNAL ARTICLE)

LA German

FS Priority Journals

EM 197709

ED Entered STN: 19900314

Last Updated on STN: 19900314

Entered Medline: 19770922

AB After phenol-acetic-acid extraction the following *Haemophilus* strains resp. their proteins were subjected the polyacrylamid-gel-electrophoresis in presence of 8 M urea: Strains of the serovar A of *H. paragallinarum*: 0083, 1516, 1598, 2213, 1645, 1646, Lohren, 2671, 1385, 758, 17756; strains of serovar B of *H. paragallinarum*: 0222, 2600, 733, 2028, 1596, 2026, 1676, 245, the S and R-form of 2403 as well as the strains 782 and 1655, which were not serotyped; strains of *H. paravium* sp. nova (HINZ: Inst. J. Syst. Bacteriol. in press): 1762, 62 (Serovar 1), 2654, 2659 (Serovar 2), 780 (Serovar 3), 94 (Serovar 4) and 1254, 0002, 0003, which were not serotyped. *H. parainfluenzae* (NCTC 4101) and *H. parasuis* were examined in the same way. The Coomassie Blue-stained protein patterns

show

that each of the strains tested developed its characteristic protein pattern, with exception of the S- and R-form of the strain 2403, which developed identical pattern. Interrelations between electrophoretic pattern and biological properties such as biochemical activities or pathogenicity could not be proved. However, the procedure described seems to be suitable for strain- or clon-identification on the subspecies

level.

The electrophoresis **apparatus**, which was made according to our

instructions was less expensive than corresponding available equipments and proved to be usable for the polyacrylamid-gel-electrophoresis.

CT *Bacterial Proteins: AN, analysis
 *Electrophoresis, Polyacrylamide Gel
 *Haemophilus: AN, analysis
 Haemophilus: CL, classification
 Species Specificity

L18 ANSWER 15 OF 15 MEDLINE
 AN 77110823 MEDLINE
 DN 77110823 PubMed ID: 827984
 TI [Qualitative and quantitative researches into bacterial flora of respiratory **apparatus** (author's transl)].
 Indagini qualitative e quantitative sulla flora batterica dell'apparato respiratorio.
 AU Torelli P C; Tortoli E
 SO ANNALI SCLAVO, (1976 Mar-Apr) 18 (2) 198-206.
 Journal code: 65C; 2985177R. ISSN: 0003-472X.
 CY Italy
 DT Journal; Article; (JOURNAL ARTICLE)
 LA Italian
 FS Priority Journals
 EM 197703
 ED Entered STN: 19900313
 Last Updated on STN: 19900313
 Entered Medline: 19770331
 AB The first group of data illustrated in this research concern the isolation's frequencies of bacterial strains from 7000 samples of respiratory materials examined in two years, 1973 and 1974, by two different techniques. Significant differences were observed in the isolation's percentage of various bacterial strains; these differences, at least partly, can be attributed at the different techniques employed in two years. The second group of data is represented by results of bacterial quantitation performed in 100 samples of purulent sputum, with respective clinical notes unknown. Potentially pathogenous strains in pure culture and at a concentration of 10(7) or higher were isolated only in 16 samples of sputum in the other samples normal bacterial flora was isolated alone or associated with potentially pathogenous strains.

CT Check Tags: Human
 Candida albicans: IP, isolation & purification
 Escherichia coli: IP, isolation & purification
Haemophilus: IP, isolation & purification
 Klebsiella: IP, isolation & purification
 Neisseria: IP, isolation & purification
 Pseudomonas aeruginosa: IP, isolation & purification
 *Respiratory System: MI, microbiology
 *Sputum: MI, microbiology
 Staphylococcus: IP, isolation & purification

=> d his

(FILE 'HOME' ENTERED AT 07:38:54 ON 09 AUG 2001)

FILE 'BIOSIS' ENTERED AT 07:39:14 ON 09 AUG 2001

L1 533 S ICT
 L2 307 S IMMUNOCHROMATOGRAPHY
 L3 31 S L1 AND L2
 L4 809 S L1 OR L2
 L5 18324 S HAEMOPHILUS OR H INFLUENZA?
 L6 1 S L5 AND L4
 L7 3855 S (POLYSACCHARIDE# OR CARBOHYDRATE#) (4A) ANTIGEN#
 L8 487656 S ANTIBODY?
 L9 1745 S L7 AND L8
 L10 0 S L9 AND L4
 L11 1 S L7 AND L4
 L12 50774 S L8 (5A) ANTIGEN#
 L13 18 S L12 AND L4
 L14 36386 S STRIP? OR TESTSTRIP
 L15 4 S L13 AND L14
 L16 ~~597 S ALL 1~~

=> d all 16;d all 111;d bib ab it 115 1-4

L6 ANSWER 1 OF 1 BIOSIS COPYRIGHT 2001 BIOSIS
 AN 2001:351661 BIOSIS
 DN PREV200100351661
 TI Laboratory diagnosis of ocular adenovirus infection: Is there really one best test.
 AU Kowalski, R. P. (1); Suzow, J. (1); Karenchak, L. M. (1); Romanowski, E. G. (1); Weck, K. E. (1); Gordon, Y. J. (1)
 CS (1) Charles T. Campbell Ophthalmic Lab and Dept of Molecular Diagnostics, Univ of Pittsburgh Med Ctr, Pittsburgh, PA USA
 SO IOVS, (March 15, 2001) Vol. 42, No. 4, pp. S579. print.
 Meeting Info.: Annual Meeting of the Association for Research in Vision and Ophthalmology Fort Lauderdale, Florida, USA April 29-May 04, 2001
 DT Conference
 LA English
 SL English
 CC Virology - Animal Host Viruses *33506
 General Biology - Symposia, Transactions and Proceedings of Conferences, Congresses, Review Annuals *00520
 Biochemical Studies - Nucleic Acids, Purines and Pyrimidines *10062
 Sense Organs, Associated Structures and Functions - Physiology and Biochemistry *20004
 Sense Organs, Associated Structures and Functions - Pathology *20006
 Physiology and Biochemistry of Bacteria *31000
 Medical and Clinical Microbiology - Virology *36006
 BC Adenoviridae 02601
 Herpesviridae 02612
 Pasteurellaceae 06703
 Chlamydiaceae 07121
 Micrococcaceae 07702
 IT Major Concepts
 Infection; Methods and Techniques; Sense Organs (Sensory Reception)

IT Parts, Structures, & Systems of Organisms
 eye: sensory system

IT Diseases
 ocular adenovirus infection: eye disease, laboratory diagnosis, viral disease

IT Chemicals & Biochemicals
 DNA: adenoviral

IT Methods & Equipment
 PCR [polymerase chain reaction]: DNA amplification, genetic method, in-situ recombinant gene expression detection, laboratory method, sensitivity, sequencing techniques, specificity; enzyme immunoassay [Adenoclone]: enzymatic method, laboratory method, sensitivity, specificity; **immunochromatography** [AdenoTest]: laboratory method, sensitivity, specificity; shell vial culture: culture method, laboratory method, sensitivity, specificity

IT Miscellaneous Descriptors
 processing time; Meeting Abstract

ORGN Super Taxa
 Adenoviridae: Animal Viruses, Viruses, Microorganisms; Animalia; Chlamydiaceae: Chlamydiales, Rickettsias and Chlamydias, Eubacteria, Bacteria, Microorganisms; Herpesviridae: Animal Viruses, Viruses, Microorganisms; Micrococcaceae: Gram-Positive Cocci, Eubacteria, Bacteria, Microorganisms; Pasteurellaceae: Facultatively Anaerobic Gram-Negative Rods, Eubacteria, Bacteria, Microorganisms

ORGN Organism Name
 Chlamydia (Chlamydiaceae); HSV [Herpes Simplex Virus] (Herpesviridae); **Haemophilus** (Pasteurellaceae); S. aureus [Staphylococcus aureus] (Micrococcaceae); St. pneumo [Staphylococcus pneumoniae] (Micrococcaceae); adenovirus (Adenoviridae): pathogen; animal (Animalia): host

ORGN Organism Superterms
 Animal Viruses; Animals; Bacteria; Eubacteria; Microorganisms; Viruses

L11 ANSWER 1 OF 1 BIOSIS COPYRIGHT 2001 BIOSIS
 AN 2001:133082 BIOSIS
 DN PREV200100133082
 TI Detection of Streptococcus pneumoniae antigen by a rapid **immunochromatographic** assay in urine samples.
 AU Dominguez, Jose (1); Gali, Nuria; Blanco, Silvia; Pedroso, Pablo; Prat, Cristina; Matas, Lurdes; Ausina, Vicente
 CS (1) Servei de Microbiologia, Hospital Universitari Germans Trias i Pujol, Carretera del Canyet s/n, 08916, Badalona, Catalonia: jadoming@ns.hugtip.scs.es Spain
 SO Chest, (January, 2001) Vol. 119, No. 1, pp. 243-249. print. ISSN: 0012-3692.
 DT Article
 LA English
 SL English
 AB Study objectives: Evaluation of a newly available rapid (15 min) **immunochromatographic** membrane test (ICT) to detect Streptococcus pneumoniae in urine samples, in order to assess its utility in the diagnosis of bacteremic and nonbacteremic pneumococcal pneumonia. Design: Retrospective study. Setting: We studied urine samples from 51 patients with bacteremic and nonbacteremic pneumonia due to S pneumoniae

diagnosed by blood culture and pneumococcal **polysaccharide** capsular **antigen** detection by counterimmunoelectrophoresis in urine samples, 16 patients with probable pneumococcal pneumonia, 71 patients with nonpneumococcal pneumonia, and 16 patients with pneumonia but no pathogen identified. Urine samples were collected and frozen at -20degreeC until used. The **ICT** test was performed following the instructions of the manufacturer. Measurements and results: S pneumoniae antigen was detected in 41 of 51 patients with pneumococcal pneumonia (80.4%); results were positive in 23 of 28 bacteremic cases (82.1%) and in 18 of 23 nonbacteremic cases (78.3%). From patients with a diagnosis of presumptive pneumococcal pneumonia, antigen was detected in seven urine samples (43.7%) and also in one case of the 16 patients with pneumonia but no pathogen identified. The specificity of the **ICT** test was 97.2%. Conclusion: The **ICT** assay is a valuable tool for the diagnosis of pneumococcal pneumonia, especially for the nonbacteremic cases.

CC Respiratory System - Pathology *16006
Clinical Biochemistry; General Methods and Applications *10006
Physiology and Biochemistry of Bacteria *31000
Medical and Clinical Microbiology - Bacteriology *36002

BC Gram-Positive Cocci 07700

IT Major Concepts
Clinical Chemistry (Allied Medical Sciences); Infection; Pulmonary Medicine (Human Medicine, Medical Sciences)

IT Diseases
non-pneumococcal pneumonia: respiratory system disease; pneumococcal pneumonia: bacterial disease, respiratory system disease

IT Chemicals & Biochemicals
Streptococcus pneumoniae antigen

IT Alternate Indexing
Pneumonia, Pneumococcal (MeSH)

IT Methods & Equipment
immunochromatographic membrane test: detection method, serodiagnostic method

ORGN Super Taxa
Gram-Positive Cocci: Eubacteria, Bacteria, Microorganisms; Hominidae: Primates, Mammalia, Vertebrata, Chordata, Animalia

ORGN Organism Name
Streptococcus pneumoniae (Gram-Positive Cocci): pathogen; human (Hominidae): patient

ORGN Organism Superterms
Animals; Bacteria; Chordates; Eubacteria; Humans; Mammals; Microorganisms; Primates; Vertebrates

L15 ANSWER 1 OF 4 BIOSIS COPYRIGHT 2001 BIOSIS
AN 2000:497615 BIOSIS
DN PREV200000497736
TI Development of rapid one-step **immunochromatographic** assay.
AU Paek, Se-Hwan (1); Lee, Seung-Hwa; Cho, Joung-Hwan; Kim, Young-Sang
CS (1) Graduate School of Biotechnology, Korea University, 5-ka, Anam-dong, Sungbuk-ku, Seoul, 136-701 South Korea
SO Methods (Orlando), (September, 2000) Vol. 22, No. 1, pp. 53-60. print.

ISSN: 1046-2023.

DT Article

LA English

SL English

AB An analytical system for a one-step immunoassay has been constructed using

the concept of **immuno chromatography**. The system employed two different antibodies that bound distinct epitopes of an analyte molecule: an antibody labeled with a signal generator (e.g., colloidal gold), which was placed in the dry state at a predetermined site on a glass-fiber membrane, and another antibody immobilized on the surface of a nitrocellulose membrane. Three membranes, one with the tracer, one with immobilized antibody, and a cellulose membrane as the absorbent of medium (in a sequence from the bottom), were attached to a plastic film and cut into **strips**. Aqueous medium containing analyte absorbed from the bottom end of the immunostrip dissolved the labeled **antibody**, and the **antigen-antibody** binding complex formed was transported into the next nitrocellulose membrane by the flow caused by capillary action. The complex subsequently reacted with the immobilized antibody, which generated a signal in proportion to the analyte concentration. The convective mass transfer of the immunoreactant to the binding partner allowed the assay to be performed with no handling of reagents. The reaction, however, was carried out under nonequilibrium conditions, which resulted in decreased sensitivity as compared with assays performed in an equilibrium mode (e.g., ELISA). To minimize such sacrifice, major factors that control system performance were identified and the system was then devised under optimal conditions.

IT Major Concepts

Methods and Techniques

IT Chemicals & Biochemicals

antibodies

IT Methods & Equipment

ELISA: analytical method, detection/labeling techniques;

immuno chromatography: analytical method, chromatographic

techniques; nitrocellulose membrane: Millipore, equipment

L15 ANSWER 2 OF 4 BIOSIS COPYRIGHT 2001 BIOSIS

AN 1998:367509 BIOSIS

DN PREV199800367509

TI Rapid accurate field diagnosis of Indian visceral leishmaniasis.

AU Sundar, Shyam; Reed, Steven G.; Singh, Vijay P.; Kumar, Prasanna C. K.; Murray, Henry W. (1)

CS (1) Box 130, 535 E. 68th St., New York, NY 10021 USA

SO Lancet (North American Edition), (Feb. 21, 1998) Vol. 351, No. 9102, pp. 563-565.

ISSN: 0099-5355.

DT Article

LA English

AB Background: A firm diagnosis of visceral leishmaniasis (kala-azar) requires demonstration of the parasite in organ aspirates or tissue

biopsy

samples. The aim of this prospective study was to assess the diagnostic usefulness of non-invasive testing for **antibody** to the leishmanial **antigen** K39 by means of antigen-impregnated nitrocellulose paper **strips** adapted for use under field conditions. Methods: One drop of peripheral blood is applied to the nitrocellulose **strip**, Three drops of test buffer

(phosphate-buffered saline plus bovine serum albumin) are added to the dried blood. The development of two visible bands indicates presence of IgG anti-K39. 323 consecutive patients with suspected kala-azar referred to two specialist units in India, and 25 healthy controls, provided fingerstick blood samples for the test. Spleen aspirates were taken from 250 patients. Findings: Kala-azar was confirmed by microscopy of spleen-aspirate smears in 127 patients. The K39 **strip** test was positive in all 127; the estimated sensitivity was therefore 100% (95%

CI

98-100). Four patients had positive **strip** tests but negative aspirate smears; all four responded to treatment for leishmaniasis. 217 individuals, including the 25 healthy controls, 73 patients with malaria or tuberculosis, and 119 spleen-aspirate-negative patients who had presumed malaria or cirrhosis (79) or no final diagnosis (40), had negative **strip**-test results. None of the 119 aspirate-negative patients developed evidence of kala-azar during 3-6 months of follow-up. The estimated specificity of the **strip** test was 98% (95-100; 217/221). Interpretation: Detection of anti-K39 by **immunochromatographic strip** testing is a rapid and non-invasive method of diagnosing kala-azar, which has good sensitivity and specificity and is well suited for use in field conditions.

IT Major Concepts

Clinical Immunology (Human Medicine, Medical Sciences); Methods and Techniques; Parasitology

IT Diseases

visceral leishmaniasis [kala-azar]: integumentary system disease, parasitic disease

IT Chemicals & Biochemicals

immunoglobulin G; recombinant K39 antigen

IT Methods & Equipment

anti-K39 antigen detection: diagnostic method, immunological method; **immunochromatographic strip** test: diagnostic method, field method, rapid, immunological method

GT India (Oriental region)

ORGN Super Taxa

Flagellata: Protozoa, Invertebrata, Animalia; Hominidae: Primates, Mammalia, Vertebrata, Chordata, Animalia

ORGN Organism Name

human (Hominidae): parasite host, patient; Leishmania-chagasi (Flagellata): parasite; Leishmania-donovani (Flagellata): parasite; Leishmania-infantum (Flagellata): parasite

ORGN Organism Superterms

Animals; Chordates; Humans; Invertebrates; Mammals; Microorganisms; Primates; Protozoans; Vertebrates

L15 ANSWER 3 OF 4 BIOSIS COPYRIGHT 2001 BIOSIS

AN 1996:285319 BIOSIS

DN PREV199699007675

TI Evaluation of **immunochromatographic** assay systems for rapid detection of hepatitis B surface **antigen** and **antibody**, Dainascreen HBsAg and Dainascreen Ausab.

AU Sato, Kumiko; Ichiyama, Satoshi (1); Iinuma, Yoshitsugu; Nada, Toshi; Shimokata, Kaoru; Nakashima, Nobuo

CS (1) Dep. Clin. Lab. Med., Nagoya Univ. Hosp., 65 Tsurumai-cho, Showa-ku, Nagoya 466 Japan

SO Journal of Clinical Microbiology, (1996) Vol. 34, No. 6, pp. 1420-1422. ISSN: 0095-1137.

DT Article
 LA English
 AB We evaluated two **immunochromatographic** assays (ICAs),
 Dainascreen HBsAg for detecting human hepatitis B surface antigen (HBsAg)
 and Dainascreen Ausab for detecting human hepatitis B surface antibody
 (anti-HBs) in human serum. The ICA systems are composed of a comb-shaped
 device that contains nitrocellulose **strips** on which complexes of
 HBsAg and anti-HBs can be visualized. The results can be read within 15
 min of incubation. The limit of detection for HBsAg was 3.1 ng/ml, and
 that for anti-HBs was 42 mIU/ml. Results of HBsAg detection agreed
 completely with those of conventional enzyme immunoassays (EIAs) and
 showed a 100% sensitivity (158 of 158 samples) and a 100% specificity
 (304 of 304 samples). The Dainascreen Ausab detected 184 of the 199
 EIA-positive samples (sensitivity, 92.5%) and yielded 6 positive results
 among the 281 EIA-negative samples (specificity, 97.9%). The ICA systems
 are rapid and sensitive methods for detecting HBsAg and anti-HBs. They
 are low-cost systems that need no complex instrumentation for analysis and
 can be recommended for routine use in clinical microbiology laboratories.

IT Major Concepts
 Immune System (Chemical Coordination and Homeostasis); Infection;
 Pathology; Serology (Allied Medical Sciences)

IT Miscellaneous Descriptors
 DIAGNOSTIC METHOD; HUMAN SERUM

ORGN Super Taxa
 Hepadnaviridae: Viruses; Hominidae: Primates, Mammalia, Vertebrata,
 Chordata, Animalia

ORGN Organism Name
 Hepadnaviridae (Hepadnaviridae); Hominidae (Hominidae)

ORGN Organism Superterms
 animals; chordates; humans; mammals; microorganisms; primates;
 vertebrates; viruses

L15 ANSWER 4 OF 4 BIOSIS COPYRIGHT 2001 BIOSIS
 AN 1985:162782 BIOSIS
 DN BR29:52778
 TI **STRIP** ASSAY FOR QUANTITATION OF **ANTIBODY** TO HEPATITIS
 B SURFACE **ANTIGEN**.
 AU VAN HAMONT J E; VINCENT J W; PAPPAS M G; HAJKOWSKI R; SETTERSTROM J A
 CS U.S. ARMY INSTITUTE OF DENTAL RESEARCH, WASHINGTON, D.C.
 SO 85TH ANNUAL MEETING OF THE AMERICAN SOCIETY FOR MICROBIOLOGY, LAS VEGAS,
 NEV., USA, MAR. 3-7, 1985. ABSTR ANNU MEET AM SOC MICROBIOL. (1985) 85
 (0), 390.
 CODEN: ASMACK. ISSN: 0094-8519.

DT Conference
 FS BR; OLD
 LA English
 IT Miscellaneous Descriptors
 ABSTRACT HUMAN MOUSE RADIOIMMUNOASSAY **IMMUNOCHROMATOGRAPHY**

=> d his

(FILE 'BIOSIS' ENTERED AT 07:39:14 ON 09 AUG 2001)
DEL HIS Y

FILE 'STNGUIDE' ENTERED AT 07:48:38 ON 09 AUG 2001

FILE 'WPIDS' ENTERED AT 07:54:13 ON 09 AUG 2001

L1 178 S IMMUNOCHROMATOGRAP? OR ICT OR IMMUNO CHROMATOGRAP?
L2 17079 S (CARBOHYDRATE? OR POLYSACCAHARIDE# OR SACCHARIDE#)
L3 28904 S L2 OR POLYSACCHARIDE#
L4 462 S L3 (5A) ANTIGEN#
L5 2 S L4 AND L1

=> d .wp 1-2

L5 ANSWER 1 OF 2 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD
AN 2000-283451 [24] WPIDS
DNC C2000-085570
TI Obtaining cell wall C-polysaccharide antigens
containing not more than 10% protein from the bacterium Streptococcus
pneumoniae for the production of antigen-specific antibodies.
DC B04 D16
IN FENT, M K; KOULCHIN, V A; MOLOKOVA, E V; MOORE, N J
PA (BINA-N) BINAX INC
CYC 36
PI WO 2000016803 A1 20000330 (200024)* EN 38p
RW: AT BE CH CY DE DK ES FI FR GB GR IE IT LU MC NL PT SE
W: AT AU CA CH CN CZ DE DK ES FI GB HU IL IN JP KR LU MX NO NZ PL PT
RU SE SK UA ZA
AU 9961513 A 20000410 (200035)
EP 1113817 A1 20010711 (200140) EN
R: AT BE CH CY DE DK ES FI FR GB GR IE IT LI LU MC NL PT SE
ADT WO 2000016803 A1 WO 1999-US21505 19990920; AU 9961513 A AU 1999-61513
19990920; EP 1113817 A1 EP 1999-948305 19990920, WO 1999-US21505 19990920
FDT AU 9961513 A Based on WO 200016803; EP 1113817 A1 Based on WO 200016803
PRAI US 1999-397110 19990916; US 1998-156486 19980918
AB WO 200016803 A UPAB: 20000522
NOVELTY - Obtaining a cell wall C-polysaccharide antigen
containing not more than 10% protein from the bacterium Streptococcus
pneumoniae for the production of antigen-specific antibodies, is new.
DETAILED DESCRIPTION - Obtaining a cell wall C-polysaccharide
antigen containing not more than 10% protein from the bacterium
Streptococcus pneumoniae comprises:
(a) culturing the bacterium for a time requisite to obtain a sample
of desired size and harvesting the bacterial cells in the form of a wet
pellet;
(b) suspending the wet cell pellet in an alkaline solution and
mixing;
(c) adjusting the pH to an acid pH with a strong acid and
centrifuging;
(d) separating a supernatant and adjusting its pH to approximate
neutrality;
(e) digesting this product with a broad spectrum protease enzyme

preparation to destroy residual proteins;

(f) adjusting the pH to the alkaline side with a weakly alkaline aqueous solution;

(g) separating out the essentially protein free **carbohydrate** or **polysaccharide antigen** on a size exclusion column equilibrated with a weakly alkaline solution; and

(h) pooling material eluted in the first peak and adjusting its pH to approximate neutrality.

INDEPENDENT CLAIMS are also included for the following:

(1) a cell wall **C-polysaccharide antigen** obtained from the method above;

(2) a method for purifying raw antibodies to *S. pneumoniae* comprising:

(a) separating from *S. pneumoniae* bacteria a cell wall **C-polysaccharide antigen** containing not more than 10% protein;

(b) conjugating the antigen to one end of a two-ended spacer molecule;

(c) coupling the conjugate to an activated chromatographic column;

(d) subjecting the raw antibodies to affinity chromatography on the column from (c) to obtain purified antigen-specific antibodies; and

(e) eluting from the column the purified antigen-specific antibodies;

(3) purified antigen-specific antibodies obtained by (2);

(4) a chromatographic column for affinity purification of raw antibodies to *S. pneumoniae* having, via a spacer molecule, a coupled purified **C-polysaccharide** cell wall **antigen** of *S. pneumoniae* containing not more than 10% protein;

(5) a method of assaying for the presence of *S. pneumoniae* or its cell wall **C-polysaccharide antigen** in a fluid comprising:

(a) extracting from a culture of *S. pneumoniae* bacteria the cell wall **C-polysaccharide antigen** containing not more than 10% protein;

(b) coupling the antigen to a spacer molecule to form a conjugate;

(c) coupling the conjugate to a chromatographic affinity column;

(d) purifying raw antibodies to *S. pneumoniae* bacteria with the chromatographic affinity column of (c) to produce purified antigen-specific antibodies; and

(e) using the purified antibodies of (d) to detect the presence/absence of *S. pneumoniae* or its **C-polysaccharide** cell wall **antigen** in a fluid;

(6) an **immuno**chromatographic (ICT) assay for the detection of *S. pneumoniae* bacteria or the **C-polysaccharide** cell wall **antigen** of the bacteria comprising:

(a) contacting a sample of a fluid suspected of containing the bacteria or their antigen with an ICT device comprising a strip of a bibulous material comprising:

(i) a first zone in which has been embedded a conjugate of a labeling

agent that displays a visible color change upon reaction of antibodies with their corresponding antigenic binding partner and purified **antigen-specific** antibodies to the **C-polysaccharide** cell wall **antigen** of *S. pneumoniae*, the antibodies having been purified by passage over a chromatographic affinity column to which is

conjugated a purified C-polysaccharide cell wall antigen of *S. pneumoniae* containing not more than 10% protein; and

(ii) a second zone having bound the same purified antigen-specific antibodies in unconjugated form, which zone is equipped with a window for viewing color changes;

the (b) allowing the sample to flow laterally along the test strip to first zone;

(c) allowing the sample, together with the conjugate of antigen-specific antibodies and label, to flow laterally along the test strip to the second zone; and

(d) within 15 minutes from the commencement of step (a) observing through the window whether a line of color has appeared in the second zone, indicating the presence in the sample of *S. pneumoniae* or its cell wall C-polysaccharide antigen; and

(7) an ICT device used in (6).

USE - The method is useful for generating antigens which, when conjugated to a spacer molecule, may be attached to a chromatographic affinity column. The column is then used to purify raw antibodies to *Streptococcus pneumoniae*. The immunochromatographic assay is used for rapid and reliable diagnosis of pathogenic states caused by *S. pneumoniae* including pneumonia, bronchitis, otitis media, sinusitis and meningitis. Dwg.0/4

L5 ANSWER 2 OF 2 WPIDS COPYRIGHT 2001 DERWENT INFORMATION LTD

AN 2000-237549 [20] WPIDS

DNC C2000-072236

TI Obtaining protein-free **carbohydrate** or **polysaccharide antigen** from a bacterium useful for immunological assays for the detection of *Legionella* caused diseases.

DC B04 D16

IN KOULCHIN, V A; MOLOKOVA, E V; MOORE, N J

PA (BINA-N) BINAX INC

CYC 36

PI WO 2000010584 A1 20000302 (200020)* EN 31p

RW: AT BE CH CY DE DK ES FI FR GB GR IE IT LU MC NL PT SE

W: AT AU CA CH CN CZ DE DK ES FI GB HU IL IN JP KR LU MX NO NZ PL PT RU SE SK UA ZA

AU 9956933 A 20000314 (200031)

EP 1107773 A1 20010620 (200135) EN

R: AT BE CH CY DE DK ES FI FR GB GR IE IT LI LU MC NL PT SE

ADT WO 2000010584 A1 WO 1999-US19506 19990825; AU 9956933 A AU 1999-56933

19990825; EP 1107773 A1 EP 1999-943941 19990825, WO 1999-US19506 19990825

FDT AU 9956933 A Based on WO 200010584; EP 1107773 A1 Based on WO 200010584

PRAI US 1998-139720 19980825

AB WO 200010584 A UPAB: 20000426

NOVELTY - A method (I) for obtaining a protein free **carbohydrate** or **polysaccharide antigen** from a known *Legionella* bacterium or serogroup of species, is new.

DETAILED DESCRIPTION - A method (I) for obtaining a protein free **carbohydrate** or **polysaccharide antigen** from a known *Legionella* bacterium or serogroup of species, is new and comprises:

(1) culturing the bacterium to obtain a sample of desired size and harvesting the bacterial cells in the form of a wet cell pellet;

(2) suspending the wet cell pellet in an alkaline solution and mixing;

(3) adjusting the pH to an acid pH with a strong acid and

centrifuging;

(4) separating the supernatant and adjusting its pH to approximate neutrality;

(5) digesting the product with a broad spectrum protease enzyme preparation to destroy residual proteins;

(6) adjusting the pH to alkaline with a weak alkaline solution; and

(7) pooling material eluted in the first peak and adjusting its pH

to

approximate neutrality.

INDEPENDENT CLAIMS are also included for the following:

(1) a protein free **carbohydrate** or **polysaccharide antigen** of a bacterium obtained from a known Legionella pneumophila species or serogroup of species using (I);

(2) a method (II) for the purification of raw antibodies to a species

or serogroup of a species of Legionella bacteria, comprising:

(a) separating from the same species or serogroup of a species of Legionella bacteria, a protein free **carbohydrate** or **polysaccharide antigen**;

(b) conjugating the antigen to one end of the a two ended spacer molecule to form a conjugate of the protein free antigen with the spacer molecule;

(c) coupling the conjugate to an activated chromatographic column;

(d) subjecting the raw antibodies to affinity chromatography on the column from (c) to obtain purified antigen specific antibodies; and

(e) eluting from the column the purified antigen-specific antibodies;

(3) purified antigen-specific antibodies to a bacterium of at least one species or serogroup of species of L. pneumophila obtained by (II);

(4) a chromatographic column (III) for affinity purification of raw antibodies in method (II);

(5) a method (IV) for assaying for the presence of Legionella bacteria or their antigenic components in a fluid comprising purifying raw

antibodies specific to L. pneumophila species and using the purified antibodies of to detect the presence or absence of the corresponding Legionella bacteria or their antigens in a fluid;

(6) a process (V) for detecting the presence of at least one species or serogroup of a species of Legionella bacteria or its antigen in a fluid

medium, where the detecting agent is antigen specific Legionella antibodies obtained by (II); and

(7) an **ICT (immunochromatographic test)** assay

(VI) for the detection of at least one species or serogroup of a species of Legionella bacteria or antigens of the bacteria, comprising:

(a) contacting a sample of a fluid suspected of containing the bacteria or their antigen with an **ICT** device comprising a strip of bibulous material comprising:

(i) a zone with a conjugate of a labelling agent (A) embedded in it, that displays a visible color change upon reaction of antibodies with their corresponding antigenic binding partner and purified antigen-specific conjugated antibodies (B) to the Legionella species to be

detected; and

(ii) a second zone having bound the same purified antigen specific antibodies in unconjugated form which is equipped with a window for viewing color changes;

the (b) allowing the sample to flow laterally along the test strip to first zone;

(c) allowing the sample together with the conjugate of affinity purified antibodies and label to flow laterally along the test strip to the second zone and;

(d) within approximately 15 minutes from the commencement of (a), observing through the window whether a line of color has appeared in the second zone thereby indicating the presence in the sample of the Legionella bacteria species or serogroup that is sought to be detected.

USE - The purified antigens are useful for the affinity purification of polyvalent antibodies to corresponding Legionella organisms. The methods are also useful for the detection of Legionella caused diseases such as Legionnaires disease and Pontiac fever in humans and for the detection of environmental sources of Legionella infectious agents.

15 ADVANTAGE - The ICT has the ability to give a test result within a minute time span for the presence or absence of L. pneumophila serogroup

1 (or its antigen) which results in high specificity and sensitivity. The test can be reliably and quickly conducted to yield a result of high specificity and high sensitivity and this is believed to be due to the strongly reactive nature of the affinity purified antibodies.
Dwg.0/2

=> d his

(FILE 'MEDLINE' ENTERED AT 07:58:03 ON 09 AUG 2001)

DEL HIS Y
E ANTIGENS, BACTERIAL+NT/CT
L1 68803 S ANTIGENS, BACTERIAL+NT/CT
L2 884 S ICT OR IMMUNOCHROMATOGRAPHY
L3 19 S L1 AND L2
E HAEMOPHILUS/CT
E E3+ALL
L4 15898 S PASTEURELLACEAE+NT/CT
L5 0 S L4 AND L3

includes all Haemophilus species

=> d .med 1-19 13

L3 ANSWER 1 OF 19 MEDLINE
AN 2001382101 MEDLINE
DN 21145371 PubMed ID: 11248521
TI On-site diagnosis of H. pylori infection by urine.
AU Miwa H; Akamatsu S; Tachikawa T; Sogabe T; Ohtaka K; Nagahara A; Sugiyama Y; Sato N
CS Department of Gastroenterology, Juntendo University School of Medicine, 2-1-1 Hongo, Bunkyo-ku, 113-8421, Tokyo, Japan..
miwahgi@med.juntendo.ac.jp
SO DIAGNOSTIC MICROBIOLOGY AND INFECTIOUS DISEASE, (2001 Feb) 39 (2) 95-7.
Journal code: DMI; 8305899. ISSN: 0732-8893.
CY United States
DT (EVALUATION STUDIES)
Journal; Article; (JOURNAL ARTICLE)
LA English
FS Priority Journals
EM 200107
ED Entered STN: 20010709
Last Updated on STN: 20010709
Entered Medline: 20010705
AB We have recently developed an on-site diagnostic kit for H. pylori infection using urine (utilizing **immunochromatographic** method employing a nitrocellulose membrane coated by extracted H. pylori antigen). Accordingly, we investigated its usefulness in 155 consecutive dyspeptic patients using the 13C urea breath test as a gold standard and further compared its performance with two commercially available rapid diagnostic kits that use whole blood (Helisal Rapid Blood, and ImmunoCard H. pylori). As the results, the urine based on-site diagnostic kit provided 95.9% sensitivity and 87.9% specificity with 92.9% accuracy, which were comparable or even better than that of both rapid whole blood tests, suggesting its usefulness in screening of H. pylori infection.
CT Check Tags: Human
Antigens, Bacterial: BL, blood
Breath Tests
*Helicobacter Infections: DI, diagnosis
*Helicobacter Infections: MI, microbiology
*Helicobacter pylori: IP, isolation & purification
Reagent Kits, Diagnostic
Sensitivity and Specificity

Urea: ME, metabolism
 *Urine: MI, microbiology

L3 ANSWER 2 OF 19 MEDLINE
 AN 2001323171 MEDLINE
 DN 21137968 PubMed ID: 11238230
 TI Early diagnosis of scrub typhus with a rapid flow assay using recombinant major outer membrane protein antigen (r56) of Orientia tsutsugamushi.
 AU Ching W M; Rowland D; Zhang Z; Bourgeois A L; Kelly D; Dasch G A; Devine P
 L
 CS Viral and Ricksettial Diseases Department, Infectious Diseases Directorate, Code 41, Naval Medical Research Center, Bethesda, Maryland, USA.. Chingw@nmrc.navy.mil
 SO CLINICAL AND DIAGNOSTIC LABORATORY IMMUNOLOGY, (2001 Mar) 8 (2) 409-14. Journal code: CB7; 9421292. ISSN: 1071-412X.
 CY United States
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 200106
 ED Entered STN: 20010611
 Last Updated on STN: 20010611
 Entered Medline: 20010607
 AB The variable 56-kDa major outer membrane protein of Orientia tsutsugamushi is the immunodominant antigen in human scrub typhus infections. We developed a rapid **immunochromatographic** flow assay (RFA) for the detection of immunoglobulin M (IgM) and IgG antibodies to O. tsutsugamushi. The RFA employs a truncated recombinant 56-kDa protein from the Karp strain as the antigen. The performance of the RFA was evaluated with a panel of 321 sera (serial bleedings of 85 individuals suspected of scrub typhus) which were collected in the Pescadore Islands, Taiwan, from 1976 to 1977. Among these 85 individuals, IgM tests were negative for 7 cases by both RFA and indirect fluorescence assay (IFA) using Karp whole-cell antigen. In 29 cases specific responses were detected by the RFA earlier than by IFA, 44 cases had the same detection time, and 5 cases were detected earlier by IFA than by RFA. For IgG responses, 4 individuals were negative with both methods, 37 cases exhibited earlier detection by RFA than IFA, 42 cases were detected at the same time, and 2 cases were detected earlier by IFA than by RFA. The sensitivities of RFA detection of antibody in sera from confirmed cases were 74 and 86% for IgM and IgG, respectively. When IgM and IgG results were combined, the sensitivity was 89%. A panel of 78 individual sera collected from patients with no evidence of scrub typhus was used to evaluate the specificity of the RFA. The specificities of the RFA were 99% for IgM and 97% for IgG. The sensitivities of IFA were 53 and 73% for IgM and IgG, respectively, and were 78% when the results of IgM and IgG were combined. The RFA test was significantly better than the IFA test for the early detection of antibody to scrub typhus in primary infections, while both tests were equally sensitive with reinfected individuals.
 CT Check Tags: Human; Support, U.S. Gov't, Non-P.H.S.

Antibody Specificity

Antigens, Bacterial: IM, immunology

*Bacterial Outer Membrane Proteins: IM, immunology

IgG: BL, blood

IgM: BL, blood

*Immunoassay: MT, methods

Immunodominant Epitopes: IM, immunology

Orientia tsutsugamushi: IM, immunology

*Orientia tsutsugamushi: IP, isolation & purification

*Reagent Strips

Recombinant Proteins: IM, immunology

*Scrub Typhus: DI, diagnosis

Sensitivity and Specificity

L3 ANSWER 3 OF 19 MEDLINE
 AN 2001318805 MEDLINE
 DN 21284474 PubMed ID: 11393288
 TI Detection of Legionella pneumophila antigen in urine samples by the
 BinaxNOW **immunochromatographic** assay and comparison with both
 Binax Legionella Urinary Enzyme Immunoassay (EIA) and Biotest Legionella
 Urin Antigen EIA.
 AU Helbig J H; Uldum S A; Luck P C; Harrison T G
 CS Medical Microbiology and Hygiene Institute, Technical University Dresden,
 Germany.. Juergen.Helbig@mailbox.tu-dresden.de
 SO JOURNAL OF MEDICAL MICROBIOLOGY, (2001 Jun) 50 (6) 509-16.
 Journal code: J2N; 0224131. ISSN: 0022-2615.
 CY England: United Kingdom
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 200106
 ED Entered STN: 20010625
 Last Updated on STN: 20010625
 Entered Medline: 20010621
 AB The new BinaxNOW **Immunochromatographic** (ICT) Assay for
 the detection of Legionella pneumophila antigens was used to test 535
 urine specimens from patients with and without Legionnaires' disease. The
 specificity, calculated by testing 112 samples from patients with
 pneumonia of aetiologies other than Legionella infection, and 167 urine
 specimens from urinary tract infections, was found to be 97.1% if the
 manufacturer's guidelines were followed. However, it was determined that
 the 'false positive' results characterised by very weak bands could be
 discounted by re-examination of the results at 60 min, yielding a
 specificity of 100%. With this minor modification of the procedure
 applied
 to examination of urine samples from 117 patients with legionellosis
 confirmed by isolation of L. pneumophila and 70 patients who had
 seroconverted to L. pneumophila serogroup 1, sensitivity was calculated
 to
 be 79.7%. In comparison, the sensitivities of the Binax Urinary Antigen
 Enzyme Immunoassay (EIA) and Biotest Urin Antigen EIA were estimated to
 be
 79.1 and 83.4%, respectively. Eleven cases (5.9%) were positive by
 BinaxNOW assay but negative by Binax or Biotest EIA, or both. The
 sensitivities of all assays increased to c. 94% if only diagnosis of
 cases
 confirmed by isolation of serogroup 1 L. pneumophila was considered,

although the sensitivity for infections caused by *L. pneumophila* serogroup

1 monoclonal antibody (MAB) subgroup Bellingham was significantly lower than for other MAB subgroups. The Biotest EIA recognised 10 (45%) of the 22 cases not caused by *L. pneumophila* serogroup 1, whereas the two Binax kits detected only three each. The ICT assay BinaxNOW can be recommended as a rapid specific test for the diagnosis of Legionnaires' diseases caused by *L. pneumophila* serogroup 1, although very weak bands should be interpreted cautiously.

CT Check Tags: Comparative Study; Human
Antibodies, Monoclonal: IM, immunology

*Antigens, Bacterial: UR, urine

Chromatography: MT, methods

False Positive Reactions

*Immunoenzyme Techniques: MT, methods

Legionella pneumophila: CL, classification

**Legionella pneumophila*: IM, immunology

Legionella pneumophila: IP, isolation & purification

*Legionnaires' Disease: DI, diagnosis

Legionnaires' Disease: UR, urine

Reagent Kits, Diagnostic

Sensitivity and Specificity

Serotyping

Time Factors

L3 ANSWER 4 OF 19 MEDLINE

AN 2001155887 MEDLINE

DN 21097861 PubMed ID: 11157611

TI Detection of *Streptococcus pneumoniae* antigen by a rapid
immunochromatographic assay in urine samples.

CM Comment in: Chest. 2001 Jan;119(1):9-11

AU Dominguez J; Gali N; Blanco S; Pedrosa P; Prat C; Matas L; Ausina V

CS Servei de Microbiologia, Hospital Universitari Germans Trias i Pujol,
Badalona, Facultat de Medicina de la Universitat Autònoma de Barcelona,
Spain.. jadoming@ns.hugtip.scs.es

SO CHEST, (2001 Jan) 119 (1) 243-9.

Journal code: D1C; 0231335. ISSN: 0012-3692.

CY United States

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Abridged Index Medicus Journals; Priority Journals

EM 200103

ED Entered STN: 20010404

Last Updated on STN: 20010404

Entered Medline: 20010322

AB STUDY OBJECTIVES: Evaluation of a newly available rapid (15 min)

immunochromatographic membrane test (ICT) to detect

Streptococcus pneumoniae in urine samples, in order to assess its utility
in the diagnosis of bacteremic and nonbacteremic pneumococcal pneumonia.

DESIGN: Retrospective study. SETTING: We studied urine samples from 51
patients with bacteremic and nonbacteremic pneumonia due to *S pneumoniae*
diagnosed by blood culture and pneumococcal polysaccharide capsular
antigen detection by counterimmunoelectrophoresis in urine samples, 16
patients with probable pneumococcal pneumonia, 71 patients with
nonpneumococcal pneumonia, and 16 patients with pneumonia but no pathogen
identified. Urine samples were collected and frozen at - 20 degrees C
until used. The ICT test was performed following the

instructions of the manufacturer. MEASUREMENTS AND RESULTS: *S. pneumoniae* antigen was detected in 41 of 51 patients with pneumococcal pneumonia (80.4%); results were positive in 23 of 28 bacteremic cases (82.1%) and in 18 of 23 nonbacteremic cases (78.3%). From patients with a diagnosis of presumptive pneumococcal pneumonia, antigen was detected in seven urine samples (43.7%) and also in one case of the 16 patients with pneumonia

but no pathogen identified. The specificity of the ICT test was 97.2%. CONCLUSION: The ICT assay is a valuable tool for the diagnosis of pneumococcal pneumonia, especially for the nonbacteremic cases.

CT Check Tags: Female; Human; Male
Adolescence
Adult
Aged
Aged, 80 and over
Bacteremia: DI, diagnosis
Bacteremia: IM, immunology
Child
Child, Preschool
Middle Age
*Pneumonia, Pneumococcal: DI, diagnosis
Pneumonia, Pneumococcal: IM, immunology
*Polysaccharides, Bacterial: UR, urine
Predictive Value of Tests

L3 ANSWER 5 OF 19 MEDLINE

AN 2000421114 MEDLINE

DN 20341579 PubMed ID: 10878074

TI Rapid diagnosis of Legionnaires' disease using an immunochromatographic assay for *Legionella pneumophila* serogroup 1 antigen in urine during an outbreak in the Netherlands.

AU Wever P C; Yzerman E P; Kuijper E J; Speelman P; Dankert J

CS Department of Medical Microbiology, Academic Medical Center, University of

Amsterdam, Amsterdam, The Netherlands.. p.c.wever@amc.uva.nl

SO JOURNAL OF CLINICAL MICROBIOLOGY, (2000 Jul) 38 (7) 2738-9.

Journal code: HSH; 7505564. ISSN: 0095-1137.

CY United States

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Priority Journals

EM 200009

ED Entered STN: 20000915

Last Updated on STN: 20000915

Entered Medline: 20000901

AB A new immunochromatographic assay for rapid qualitative detection of *Legionella pneumophila* serogroup 1 antigen in urine specimens

was used during an outbreak of legionellosis in The Netherlands. The assay

seems of the utmost value in providing a rapid diagnosis of Legionnaires' disease in patients with severe community-acquired pneumonia in an outbreak setting.

CT Check Tags: Human

*Antigens, Bacterial: UR, urine

Chromatography: MT, methods
 *Disease Outbreaks
 *Immunologic Techniques
 Legionella pneumophila: IM, immunology
 *Legionella pneumophila: IP, isolation & purification
 *Legionnaires' Disease: DI, diagnosis
 Legionnaires' Disease: EP, epidemiology
 Legionnaires' Disease: MI, microbiology
 Netherlands: EP, epidemiology

L3 ANSWER 6 OF 19 MEDLINE
 AN 2000153062 MEDLINE
 DN 20153062 PubMed ID: 10691203
 TI Evaluation of a rapid **immunochromatographic** assay for the
 detection of Legionella antigen in urine samples.
 AU Dominguez J; Gali N; Matas L; Pedroso P; Hernandez A; Padilla E; Ausina V
 CS Servicio de Microbiologia, Hospital Universitari Germans Trias i Pujol,
 Facultad de Medicina de la Universitat Autonoma de Barcelona, Spain..
 jadoming@ns.hugtip.scs.es
 SO EUROPEAN JOURNAL OF CLINICAL MICROBIOLOGY AND INFECTIOUS DISEASES, (1999
 Dec) 18 (12) 896-8.
 Journal code: EM5; 8804297. ISSN: 0934-9723.
 CY GERMANY: Germany, Federal Republic of
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 200003
 ED Entered STN: 20000327
 Last Updated on STN: 20000327
 Entered Medline: 20000314
 AB A new **immunochromatographic** membrane assay for detecting
 Legionella pneumophila serogroup 1 antigen in urine samples (Binax Now
 Legionella Urinary Antigen Test; Binax, USA) was evaluated. Its
 sensitivity, specificity and level of agreement with the Binax enzyme
 immunoassay were compared using nonconcentrated and concentrated urine
 samples. The overall agreement between the two tests was 98.1%; the
 specificity of both was 100%. The sensitivity of the
 immunochromatographic assay was 55.5% in nonconcentrated urine and
 97.2% in concentrated urine in comparison with the enzyme immunoassay,
 using concentrated urine as the reference test. This
 immunochromatographic assay screens successfully for Legionella
 pneumophila serogroup 1 soluble antigen in concentrated urine samples.
 CT Check Tags: Female; Human; Male
 Adult
 Aged
 *Antigens, Bacterial: UR, urine
 *Chromatography: MT, methods
 Evaluation Studies
 Immunoenzyme Techniques
 *Immunologic Techniques
 Legionella pneumophila: IM, immunology
 *Legionella pneumophila: IP, isolation & purification
 *Legionnaires' Disease: DI, diagnosis
 Legionnaires' Disease: MI, microbiology
 Middle Age
 Sensitivity and Specificity
 *Urine: MI, microbiology

L3 ANSWER 7 OF 19 MEDLINE
 AN 2000087388 MEDLINE
 DN 20087388 PubMed ID: 10618283
 TI Detection of Francisella tularensis in biological specimens using a capture enzyme-linked immunosorbent assay, an **immunochromatographic** handheld assay, and a PCR.
 AU Grunow R; Splettstoesser W; McDonald S; Otterbein C; O'Brien T; Morgan C; Aldrich J; Hofer E; Finke E J; Meyer H
 CS Institute of Microbiology, Federal Armed Forces Medical Academy, 80937 Munich, Germany.. tb101cn@mail.lrz-muenchen.de
 SO CLINICAL AND DIAGNOSTIC LABORATORY IMMUNOLOGY, (2000 Jan) 7 (1) 86-90.
 Journal code: CB7; 9421292. ISSN: 1071-412X.
 CY United States
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 200003
 ED Entered STN: 20000314
 Last Updated on STN: 20000314
 Entered Medline: 20000302
 AB The early detection of Francisella tularensis, the causative agent of tularemia, is important for adequate treatment by antibiotics and the outcome of the disease. Here we describe a new capture enzyme-linked immunosorbent assay (cELISA) based on monoclonal antibodies specific for lipopolysaccharide (LPS) of Francisella tularensis subsp. holarctica and Francisella tularensis subsp. tularensis. No cross-reactivity with Francisella tularensis subsp. novicida, Francisella philomiragia, and a panel of other possibly related bacteria, including Brucella spp., Yersinia spp., Escherichia coli, and Burkholderia spp., was observed. The detection limit of the assay was 10(3) to 10(4) bacteria/ml. This sensitivity was achieved by solubilization of the LPS prior to the cELISA.
 In addition, a novel **immunochromatographic** membrane-based handheld assay (HHA) and a PCR, targeting sequences of the 17-kDa protein (TUL4) gene of F. tularensis, were used in this study. Compared to the cELISA, the sensitivity of the HHA was about 100 times lower and that of the PCR was about 10 times higher. All three techniques were successfully applied to detect F. tularensis in tissue samples of European brown hares (Lepus europaeus). Whereas all infected samples were recognized by the cELISA, those with relatively low bacterial load were partially or not detected by PCR and HHA, probably due to inhibitors or lack of sensitivity. In conclusion, the HHA can be used as a very fast and simple approach to perform field diagnosis to obtain a first hint of an infection with F. tularensis, especially in emergent situations. In any suspect case, the diagnosis should be confirmed by more sensitive techniques, such as the cELISA and PCR.
 CT Check Tags: Comparative Study
 Antibodies, Monoclonal
 Blotting, Western
 Brucella: IM, immunology
 Burkholderia: IM, immunology
 Cross Reactions
 *Enzyme-Linked Immunosorbent Assay
 Enzyme-Linked Immunosorbent Assay: MT, methods

Escherichia coli: IM, immunology
 Francisella tularensis: CL, classification
 *Francisella tularensis: GE, genetics
 *Francisella tularensis: IP, isolation & purification
 *Immunologic Tests: MT, methods
 *Lipopolysaccharides: IM, immunology
 *Polymerase Chain Reaction
 Polymerase Chain Reaction: MT, methods
 Sensitivity and Specificity
 Yersinia: IM, immunology

L3 ANSWER 8 OF 19 MEDLINE
 AN 1999308407 MEDLINE
 DN 99308407 PubMed ID: 10383257
 TI Rapid **immunochromatographic** assay for diagnosis of tuberculosis:
 antibodies detected may not be specific.
 CM Comment on: J Clin Microbiol. 1998 Nov;36(11):3443
 AU Freeman R; Magee J; Barratt A; Wheeler J; Steward M; Lee M; Piggott N
 SO JOURNAL OF CLINICAL MICROBIOLOGY, (1999 Jun) 37 (6) 2111-2.
 Journal code: HSH; 7505564. ISSN: 0095-1137.
 CY United States
 DT Commentary
 Letter
 LA English
 FS Priority Journals
 EM 199906
 ED Entered STN: 19990712
 Last Updated on STN: 20000303
 Entered Medline: 19990623
 CT Check Tags: Human
 Antibodies, Monoclonal
Antigens, Bacterial: IM, immunology
 Blotting, Western
 Molecular Weight
 Recombinant Proteins: IM, immunology
 *Tuberculosis: DI, diagnosis
 Tuberculosis: IM, immunology
 Tuberculosis, Pulmonary: DI, diagnosis
 Tuberculosis, Pulmonary: IM, immunology

L3 ANSWER 9 OF 19 MEDLINE
 AN 1998453872 MEDLINE
 DN 98453872 PubMed ID: 9780587
 TI Evaluation of **immunochromatography**-based rapid detection kit for
 fecal Escherichia coli O157.
 AU Takeda T; Yamagata K; Yoshida Y; Yoshino K; Nomura T
 CS Department of Infectious Diseases Research, National Children's Medical
 Research Center.
 SO KANSENSHOGAKU ZASSHI. JOURNAL OF THE JAPANESE ASSOCIATION FOR INFECTIOUS
 DISEASES, (1998 Aug) 72 (8) 834-9.
 Journal code: IJR; 0236671. ISSN: 0387-5911.
 CY Japan
 DT Journal; Article; (JOURNAL ARTICLE)
 LA Japanese
 FS Priority Journals
 EM 199812
 ED Entered STN: 19990115

Last Updated on STN: 19990115

Entered Medline: 19981201

AB "Quix" is an **immunochromatography**-based direct detection kit for the E. coli O157 LPS antigen in the patient's stool. The present study

was conducted to evaluate the efficacy of the kit for rapid diagnosis of enterohemorrhagic Escherichia coli (EHEC) O157 infection. Sensitivity of the kit was determined using a pure culture of a clinical isolate of E. coli O157. Analytical sensitivity was found to be 5×10^5 CFU/ml. When compared with the culture method using fecal samples of 64 patients and with bloody diarrhea, sensitivity and specificity were 95.0% (19/20) and 86.4% (38/44), respectively, and overall agreement to culture method was 89.1% (57/64). One patient was found positive by culture method while negative in the present method, where the sample contained a low number

of the cells less than the detection limit. Four of the six patients with a negative result by culture method and positive in the present method,

were confirmed E. coli O157 infection by positive IgM antibody response

against the E. coli O157 LPS. The discrepancy between the two methods seemed to be

attributable to antibiotic administration. In one patient, Salmonella urbana (O30(1)30(2)) was detected. The O30(1) antigen of this bacterium is

well known to be identical to the E. coli O157 antigen. When the present method was compared with an ELISA-based E. coli O157 LPS antigen detection

kit, sensitivity and specificity were 100% (11/11) and 82.1% (23/28), and overall agreement to ELISA method was 87.2% (34/39). From these findings, Quix is useful as a rapid diagnostic kit in the primary clinics, outpatient or bedside use. E. coli O157 LPS antigen in patient's fecal samples can be detected in about five minutes with this simple procedure. Early diagnosis using such a simple kit will largely contribute for the early treatment and prevention of severe complications of the E. coli

O157 infection.

CT Check Tags: Female; Human; Male

Antigens, Bacterial: AN, analysis

Child

Child, Preschool

Escherichia coli Infections: DI, diagnosis

Escherichia coli O157: IM, immunology

*Escherichia coli O157: IP, isolation & purification

Evaluation Studies

*Feces: MI, microbiology

Reagent Kits, Diagnostic

Sensitivity and Specificity

L3 ANSWER 10 OF 19 MEDLINE

AN 1998193869 MEDLINE

DN 98193869 PubMed ID: 9532690

TI [A comparative study of the information value of the immunofluorescence and **immunochromatographic** identification of Chlamydia trachomatis antigens in smear material from the cervical canal and of the cytological picture of a vaginal discharge in pregnant women].
Sravnitel'noe izuchenie informativnosti immunofluorescentnoi i

immunokhromatograficheskoi identifikatsii antigenov Chlamydia trachomatis v sosskobnom materiale iz tservikal'nogo kanala i tsitologicheskoi kartiny vaginal'nogo otdeliaemogo u beremennykh.

AU Zul'karneev R Sh; Kalinin Iu T; Afanas'ev S S; Rubal'skii O V; Denisov L A; Vorob'ev A A; Kiselev V I; Afanas'ev D S; Voronin M V; Afanas'ev M S

CS Health Department of the City Administration of Astrakhan, Russia.

SO ZHURNAL MIKROBIOLOGII, EPIDEMIOLOGII I IMMUNOBIOLOGII, (1998 Jan-Feb) (1) 64-7.

Journal code: Y90; 0415217. ISSN: 0372-9311.

CY RUSSIA: Russian Federation

DT Journal; Article; (JOURNAL ARTICLE)

LA Russian

FS Priority Journals

EM 199806

ED Entered STN: 19980618
Last Updated on STN: 19980618
Entered Medline: 19980608

AB The cytological picture of the vaginal discharge and scrape material, obtained from 30 pregnant women of the control group (not infected with C.trachomatis) and 61 pregnant women with chlamydiosis; of these, in 42 women the comparative identification of chlamydiae by the methods of direct immunofluorescence and **immunochromatography** was carried out. Direct immunofluorescence was carried out with the use of a set of reagents "MicroTrak" (USA) and **immunochromatographic** identification was made with the use of a set of reagents "Clearview Chlamydia" (Britain). The comparison of the results of **immunochromatography** and direct immunofluorescence revealed that the sensitivity of the set "Clearview Chlamydia" was 100.0% in comparison with the data obtained in the examination of women with the use of the set "MicroTrak". The negative results coincided in 90.1% of cases. The common features of the cytological picture of vaginal samples taken from pregnant women were established. The optimum system of the examination of pregnant women suspected for chlamydiosis, as well as for the evaluation of the effectiveness of its treatment, was proposed.

CT Check Tags: Comparative Study; Female; Human
*Antigens, Bacterial: AN, analysis
*Chlamydia Infections: DI, diagnosis
*Chlamydia trachomatis: IM, immunology
Chromatography
Fluorescent Antibody Technique, Direct
Immunologic Tests
Pregnancy
*Pregnancy Complications, Infectious: DI, diagnosis
Sensitivity and Specificity
Time Factors
*Vaginal Discharge: DI, diagnosis
*Vaginal Smears: MT, methods

L3 ANSWER 11 OF 19 MEDLINE

AN 1998134302 MEDLINE

DN 98134302 PubMed ID: 9474027

TI Development of FlexSure HP--an **immunochromatographic** method to detect antibodies against Helicobacter pylori.

AU Schrier W H; Schoengold R J; Baker J T; Norell J L; Jaseph C L; Okin Y; Doe J Y; Chandler H

CS SmithKline Diagnostics, Inc., Palo Alto, CA 94303, USA.
 SO CLINICAL CHEMISTRY, (1998 Feb) 44 (2) 293-8.
 Journal code: DBZ; 9421549. ISSN: 0009-9147.
 CY United States
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 199803
 ED Entered STN: 19980326
 Last Updated on STN: 19980326
 Entered Medline: 19980316
 AB We describe a solid-phase **immunochromatographic** serologic test, FlexSure HP, to detect IgG antibodies against Helicobacter pylori. H. pylori colonize the stomach and proximal duodenum, cause ulcer disease and mucosa-associated lymphoid tissue lymphoma, and have a role in the development of other disorders, including gastric adenocarcinoma.
 FlexSure HP consists of a test strip, conjugate pad, and absorbent pad, in a novel reverse-flow chromatography format. In these studies, FlexSure HP was demonstrated to be specific for IgG antibodies against H. pylori. The reactive cutoff of the test was consistent with [13C]urea breath test and commercially available ELISAs. FlexSure HP had 94% sensitivity, 88% specificity, and 91% accuracy relative to [13C]urea breath test; and 95% sensitivity, 94% specificity, and 95% overall agreement relative to high-molecular-mass cell-associated protein enzyme immunoassay (HM-CAP EIA). FlexSure HP is a simple-to-perform, visually read test requiring no specialized training, equipment, or instrumentation, and yields rapid, accurate, qualitative results.
 CT Check Tags: Female; Human; Male
 Adult
 *Antibodies: BL, blood
 Antigens, Bacterial: BL, blood
 Antigens, Bacterial: IM, immunology
 Breath Tests
 Cross Reactions: IM, immunology
 *Helicobacter Infections: BL, blood
 Helicobacter Infections: IM, immunology
 *Helicobacter pylori: IM, immunology
 Immunosorbents: ME, metabolism
 Middle Age
 ROC Curve
 Reagent Kits, Diagnostic
 Reproducibility of Results
 Sensitivity and Specificity
 Stomach: MI, microbiology
 Urea: AN, analysis
 L3 ANSWER 12 OF 19 MEDLINE
 AN 96388854 MEDLINE
 DN 96388854 PubMed ID: 8796254
 TI Clinical evaluation of a rapid **immunochromatographic** assay based on the 38 kDa antigen of Mycobacterium tuberculosis on patients with pulmonary tuberculosis in China.
 CM Comment in: Tuber Lung Dis. 1997;78(1):85
 AU Cole R A; Lu H M; Shi Y Z; Wang J; De-Hua T; Zhou A T
 CS ICT Diagnostics, Sydney, Australia.

SO TUBERCLE AND LUNG DISEASE, (1996 Aug) 77 (4) 363-8.
Journal code: A8C; 9212467. ISSN: 0962-8479.

CY SCOTLAND: United Kingdom

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Priority Journals

EM 199610

ED Entered STN: 19961106
Last Updated on STN: 19990129
Entered Medline: 19961022

AB SETTING: A rapid membrane-based antibody assay capable of diagnosing pulmonary tuberculosis within 15 min has been developed using the 38 kDa antigen from Mycobacterium tuberculosis. OBJECTIVE: To determine the specificity and sensitivity of this assay and evaluate its usefulness in

a clinical setting. DESIGN: Sera from patients with active pulmonary tuberculosis were obtained from three hospitals in China. The control groups consisted of patients who were diagnosed with lung diseases other than tuberculosis and healthy subjects. RESULTS: Antibody was detected in 54 of 61 (89%) sputum positive patients and 67 of 91 (74%) sputum negative patients who had been clinically diagnosed as having active pulmonary tuberculosis. Five out of 56 (9%) patients with respiratory diseases other than tuberculosis and 1 out of 30 (3%) healthy controls had a positive antibody response. The overall specificity of the assay was 93% and the positive predictive value was 95%. We conclude that this assay is rapid, sensitive and specific and will be a valuable aid in the clinical diagnosis of pulmonary tuberculosis.

CT Check Tags: Female; Human; Male
Adult
***Antigens, Bacterial: AN, analysis**
*Immunoassay: MT, methods
Middle Age
*Mycobacterium tuberculosis: IM, immunology
Mycobacterium tuberculosis: IP, isolation & purification
Predictive Value of Tests
Respiratory Tract Diseases: IM, immunology
Sensitivity and Specificity
Sputum: MI, microbiology
Time Factors
*Tuberculosis, Pulmonary: DI, diagnosis
Tuberculosis, Pulmonary: MI, microbiology

L3 ANSWER 13 OF 19 MEDLINE

AN 94154898 MEDLINE

DN 94154898 PubMed ID: 8111517

TI Comparison of an **immunochromatographic** method for rapid identification of group A streptococcal antigen with culture method.

AU Ehrlich T P; Schwartz R H; Wientzen R; Thorne M M

CS Department of Family Practice, Medical College of Virginia, Fairfax Family Practice Center Inc.

SO ARCHIVES OF FAMILY MEDICINE, (1993 Aug) 2 (8) 866-9.
Journal code: BX6; 9300357. ISSN: 1063-3987.

CY United States

DT Journal; Article; (JOURNAL ARTICLE)

LA English
FS Priority Journals
EM 199403

ED Entered STN: 19940406
Last Updated on STN: 19940406
Entered Medline: 19940330

AB OBJECTIVES: To compare the sensitivity and specificity of Concise Strep A (Hybritech, San Diego, Calif), an **immunochromatographic** group A streptococcal rapid antigen detection system, with a two-plate culture method for the diagnosis of streptococcal pharyngitis, and to evaluate

the need for routine back-up culture when this rapid test is used. DESIGN: Throat cultures were obtained from 351 children with acute pharyngitis by duplicate rayon-tipped swabs held in parallel and vigorously rubbed against both tonsils and the posterior pharyngeal wall. One swab was tested for group A streptococcal antigen by a registered licensed laboratory technologist in the pediatrician's office. The other swab was streaked over each of two sheep blood agar plates, one of which was enhanced with trimethoprim in combination with sulfamethoxazole. The

plain sheep blood agar plate was then incubated in a candle-extinguish jar. The enhanced agar plate was placed in a gas-pack anaerobic jar. Both plates were incubated for up to 48 hours at 35 degrees C. SETTING: A six-person group pediatric practice. PARTICIPANTS: Three hundred fifty-one children. RESULTS: The Concise Strep A antigen detection test produced 129 positive results. Only six of the 129 were not confirmed by culture method. There were four false-negative rapid streptococcal antigen detection test results, all of which were found after a single overnight incubation. The sensitivity for the Concise Strep A test was 96.9% and the specificity

was 97.4%. The plain 5% sheep blood agar plate (without trimethoprim and sulfamethoxazole), which was incubated in a candle-extinguish jar, identified 123 (97%) of the 127 positive throat cultures. The second 24-hour incubation and use of trimethoprim and sulfamethoxazole agar were not rewarding for this study. CONCLUSIONS: Concise Strep A, a polyclonal antibody test, in conjunction with a color **immunochromatographic** assay for soluble streptococcal carbohydrate antigen A appears to be accurate, sensitive, and specific when throat swabs are carefully

obtained and when qualified, licensed laboratory technologists perform the procedure. Further studies should be done to confirm our findings, especially when nurses or office staff perform the rapid test procedure

in the office setting. If our findings are confirmed, the use of back-up cultures for negative rapid test results obtained using Concise Strep A would be unnecessary.

CT Check Tags: Comparative Study; Human; Support, Non-U.S. Gov't
Acute Disease

*Antigens, Bacterial: AN, analysis

Bacteriological Techniques

Child

Chromatography

False Negative Reactions

Immunoassay

*Pharyngitis: DI, diagnosis

Pharyngitis: MI, microbiology

Pharynx: MI, microbiology

Reagent Kits, Diagnostic
Sensitivity and Specificity
*Streptococcal Infections: DI, diagnosis
Streptococcus pyogenes: IM, immunology
*Streptococcus pyogenes: IP, isolation & purification

L3 ANSWER 14 OF 19 MEDLINE
AN 93036353 MEDLINE
DN 93036353 PubMed ID: 1416045
TI Kinetic chromatographic sequential addition immunoassays using protein A affinity chromatography.
AU Cassidy S A; Janis L J; Regnier F E
CS Department of Chemistry, Purdue University, West Lafayette, Indiana 47907.
NC GM 25431 (NIGMS)
SO ANALYTICAL CHEMISTRY, (1992 Sep 1) 64 (17) 1973-7.
Journal code: 4NR; 0370536. ISSN: 0003-2700.
CY United States
DT Journal; Article; (JOURNAL ARTICLE)
LA English
FS Priority Journals
EM 199210
ED Entered STN: 19930122
Last Updated on STN: 19930122
Entered Medline: 19921026
AB A new type of chromatographic immunoassay based on sequential addition is described. On a protein A column, the antibody, the sample containing the antigen, and then a known amount of antigen are sequentially injected. This assay is designed to shorten analysis times and reduce complexity of dual-column chromatographic immunoassays, circumvent desorption buffer interferences common to affinity chromatography, and eliminate the need for tagged molecules. This new technique is named kinetic **immunochromatography** sequential addition (KICQA). Because of its kinetic nature, flow rate will have a large effect on KICQA, and the impact of changing flow rate is studied extensively. By use of various amounts of antibody, the dynamic range of KICQA is shown to be selectable over 2.5 orders of magnitude. Finally, KICQA was used to determine transferrin and albumin in human serum. Both analytes show good agreement with their respective reference methods, and an albumin assay was performed in under 1 min.
CT Check Tags: Support, U.S. Gov't, P.H.S.
Antigen-Antibody Complex
*Apoproteins: AN, analysis
*Chromatography, Affinity
*Immunoassay
*Serum Albumin: AN, analysis
***Staphylococcal Protein A: CH, chemistry**
*Transferrin: AN, analysis

L3 ANSWER 15 OF 19 MEDLINE
AN 90028644 MEDLINE
DN 90028644 PubMed ID: 2679905
TI [Monoclonal antibodies cross-reacting with fibroblasts of interstitial connective tissue of the myocardium and cell wall protein antigens of group A Streptococcus].
Monoklonal'nye antitela, perekrestno reagiruiushchie s fibroblastami interstitsial'noi soedinitel'noi tkani miokarda i belkovymi antigenami

kletochnoi stenki streptokokka gruppy A.
 AU Abyzov V N; Drobyshevskaja E I; Liampert I M; Borodiiuk N A; Panasiuk A F
 SO BIULLETEN EKSPERIMENTALNOI BIOLOGII I MEDITSINY, (1989 Jul) 108 (7) 74-6.
 Journal code: A74; 0370627. ISSN: 0365-9615.
 CY USSR
 DT Journal; Article; (JOURNAL ARTICLE)
 LA Russian
 FS Priority Journals
 EM 198912
 ED Entered STN: 19900328
 Last Updated on STN: 19900328
 Entered Medline: 19891201
 AB Monoclonal antibodies (MCA) B6/5 and C5/3 were obtained after
 immunization
 of BALB/c mice with the protein non-type-specific antigens (NTSA) of
 streptococcal group A cell wall. MCA B6/5 in the indirect
 immunofluorescence react with human and animal interstitial connective
 tissue (ICT) of the myocardium and human fibroblast culture
 cells. MCA C5/3 react with the bands of muscle fibers of the myocardium.
 MCA B6/5 and C5/3 are autoantibodies. It was revealed that these MCA are
 directed to two streptococcal cross-reacting antigens (CRA). Production
 of B6/5 and C5/3, apparently, does not depend on the possibility of some
 streptococcal antigens to bind fibrinogen. Bound immunoglobulins were not
 revealed in the ICT and in the muscle fibres by the cultivation
 of the C5/3 monoclonal. Firstly it was stated that, MCA B6/5, reacting
 with fibroblasts and with streptococcal CRA, are capable to fix in the
 ICT of myocardium, what is typical for the phenomenon described in
 rheumatic fever.
 CT Check Tags: Animal; Human
 *Antibodies, Monoclonal: IM, immunology
 *Antigens, Bacterial: IM, immunology
 *Autoantibodies: AN, analysis
 Cross Reactions
 Fibroblasts: IM, immunology
 Fluorescent Antibody Technique
 Mice
 Mice, Inbred BALB C
 *Myocardium: IM, immunology
 Rheumatic Fever: IM, immunology
 *Streptococcus pyogenes: IM, immunology
 L3 ANSWER 16 OF 19 MEDLINE
 AN 81035609 MEDLINE
 DN 81035609 PubMed ID: 6158813
 TI [T-proteins of Streptococcus pyogenes. II. Communication: Preparation of
 specific immunoabsorbents for isolation of anti-T-antibodies and
 investigations of the T4/24...pattern (author's transl)].
 T-Proteine des Streptococcus pyogenes. II. Mitteilung: Preparation
 spezifischer Immunoabsorbentien zur Isolierung von T-Antikörpern und
 Untersuchungen de T4/24...Agglutinationsgruppe.
 AU Schmidt K H; Kuhnemund O; Kohler W
 SO ZENTRALBLATT FÜR BAKTERIOLOGIE. 1. ABT. ORIGINALE. A: MEDIZINISCHE
 MIKROBIOLOGIE, INFektionsKRANKHEITEN UND PARASITIOLOGIE, (1980) 246 (4)
 489-98.
 Journal code: Y5P; 8005748. ISSN: 0172-5599.

CY GERMANY, WEST: Germany, Federal Republic of
 DT Journal; Article; (JOURNAL ARTICLE)
 LA German
 FS Priority Journals
 EM 198012
 ED Entered STN: 19900316
 Last Updated on STN: 19900316
 Entered Medline: 19801218
 AB The T4-antigen(s) in Streptococcus pyogenes standard type strains of the T4/24...pattern was (were) investigated in regard to the uniformity or diversity of the antigen(s) in this T-complex. Tryptic T-extracts of types 4, 24, 26, 28, 29, 46, 48, 60 were purified by ion exchange chromatography on DEAE-cellulose. Anti-T4-antibodies were isolated by immunoadsorption chromatography on AH-sepharose linked T4-antigen. Purified T4-antigen showed in SDS-electrophoresis a similar multiple molecular size structure as T1-antigen described earlier. Comparative serological studies of T-antigens of types 4, 24, 29 and 46 revealed reactions of identity to anti-T4-antibodies in Ouchterlony tests. Extracts of types 26, 28, 48 and 60 did not precipitate with anti-T4-antibodies, but types 28 and 48 showed crossreaction to the relevant antisera (anti-T28 and anti-T48, resp.) obviously caused by traces of R-28 antigen in both antigen preparations. Strains of the types 4, 24, 29, 46, 48 and 60 were agglutinated by anti-T4-antibodies. The reaction could be inhibited by T4-antigen. The strains of type 26 and 28 used in our experiments did not contain T4-antigen. Agglutination as well as immunoprecipitation reactions with specific antibodies prepared by **immunochromatography** proved the existence of common T4-antigenic determinants in type 4, 24, 29, 46, 48 and 60.
 CT Check Tags: Animal
 Antibodies, Bacterial: IM, immunology
 *Antibodies, Bacterial: IP, isolation & purification
 *Antigen-Antibody Complex
 *Antigens, Bacterial: IM, immunology
 Antigens, Bacterial: IP, isolation & purification
 *Bacterial Proteins: IM, immunology
 Epitopes
 Immunodiffusion
 *Immunosorbents
 Rabbits
 *Streptococcus pyogenes: IM, immunology
 L3 ANSWER 17 OF 19 MEDLINE
 AN 80243108 MEDLINE
 DN 80243108 PubMed ID: 6994834
 TI [Study of cross reactions between group A streptococcal antigens and cardiac interstitial connective tissue fibroblasts of different species of mammals].
 Izuchenie perekrestnykh reaktsii mezhdu antigenami streptokokka gruppy A i fibroblastami interstitsial'noi soedinitel'noi tkani serdtsa mlekopitaiushchikh raznykh vidov.
 AU Kochetkova E V; Liampert I M; Kolesnikova V Iu; Semenova E N
 SO BIULLETEN EKSPERIMENTALNOI BIOLOGII I MEDITSINY, (1980 May) 89 (5) 582-4.

Journal code: A74; 0370627. ISSN: 0365-9615.

CY USSR

DT Journal; Article; (JOURNAL ARTICLE)

LA Russian

FS Priority Journals

EM 198010

ED Entered STN: 19900315
Last Updated on STN: 19900315
Entered Medline: 19801027

AB It has been revealed by indirect immunofluorescence that sera of rabbits immunized with group A streptococcal non-type-specific antigens and F (ab')₂ fragments of IgG obtained from these sera intensively react with fibroblasts of interstitial connective tissue (ICT) of the myocardium of humans and different mammalian species with the exception of rabbits. Negative or weak reactions with the same sera and F (ab')₂ fragments were observed with fibroblasts of ICT of nonimmunized rabbits. Myocardial ICT of immunized animals showed bound immunoglobulins. This model can be used for making clear the reason for presence in sera of rheumatic patients circulating antibodies only to heterologous but not to homologous ICT of the myocardium regardless of the revealing bound immunoglobulins in the ICT of these patients.

CT Check Tags: Animal; Comparative Study; Human
Antigens, Bacterial: IM, immunology
Autoantibodies: AN, analysis
Cattle
*Cross Reactions
Fibroblasts: IM, immunology
Fluorescent Antibody Technique
Guinea Pigs
Heart: EM, embryology
Immunoglobulins, Fab: IM, immunology
Myocardium: CY, cytology
*Myocardium: IM, immunology
Rabbits: IM, immunology
Rats
Rheumatic Fever: IM, immunology
Species Specificity
*Streptococcus pyogenes: IM, immunology
Swine: IM, immunology

L3 ANSWER 18 OF 19 MEDLINE

AN 80181173 MEDLINE

DN 80181173 PubMed ID: 6989752

TI Mitogenic and antigenic properties of group A streptococcal M protein preparations.

AU Knoll H; Kuhnemund O; Havlicek J

SO IMMUNOBIOLOGY, (1980 Jan) 156 (4-5) 537-48.
Journal code: GH3; 8002742. ISSN: 0171-2985.

CY GERMANY, WEST: Germany, Federal Republic of

DT Journal; Article; (JOURNAL ARTICLE)

LA English

FS Priority Journals

EM 198007

ED Entered STN: 19900315
Last Updated on STN: 19900315

Entered Medline: 19800728

AB The separation of a mitogenic substance in M protein preparations of Streptococcus pyogenes, type 1 and type 12, is described. The isolation was achieved by gel chromatography on Biogel A 0.5 m, and by **immunochromatography** on immobilized type specific antibodies. In the delayed cutaneous hypersensitivity test the **immunochromatography** purified M 1 protein caused erythema but no unspecific mitogenicity could be found by lymphocyte transformation test. In neutralization experiments the mitogenic activity of M protein preparations was specifically inhibited by anti-erythrogenic toxin antisera.

CT Check Tags: Animal; Human
 Adult
 Antibody Specificity
***Antigens, Bacterial**
 *Bacterial Proteins: IM, immunology
 Bacterial Proteins: IP, isolation & purification
 Chromatography, Gel
 Fractionation
 Hydrochloric Acid: IM, immunology
 Hypersensitivity, Delayed: IM, immunology
 Middle Age
 *Mitogens: PD, pharmacology
 Neutralization Tests
 Precipitation
 Rabbits
 *Streptococcus pyogenes: IM, immunology

L3 ANSWER 19 OF 19 MEDLINE
 AN 80061440 MEDLINE
 DN 80061440 PubMed ID: 388938
 TI Immunolectron microscopic localization of T proteins in the cell wall of Streptococcus pyogenes.
 AU Wagner B; Schmidt K H; Wagner M
 SO ZENTRALBLATT FUR BAKTERIOLOGIE, PARASITENKUNDE, INFEKTIONSKRANKHEITEN UND HYGIENE. ERSTE ABTEILUNG ORIGINALE. REIHE A: MEDIZINISCHE MIKROBIOLOGIE UND PARASITOLOGIE, (1979 Jul) 244 (2-3) 192-201.
 Journal code: Y52; 0331570. ISSN: 0300-9688.
 CY GERMANY, WEST: Germany, Federal Republic of
 DT Journal; Article; (JOURNAL ARTICLE)
 LA English
 FS Priority Journals
 EM 198001
 ED Entered STN: 19900315
 Last Updated on STN: 19900315
 Entered Medline: 19800128

AB T proteins of Streptococcus pyogenes (group A streptococci) were localized by means of **immunochromatographical** isolated anti-T-antibodies. For the electron microscopical detection both the direct and the indirect immunoferritin techniques were used. The arrangement of the ferritin particles showed, that the T proteins are evenly distributed on the whole cell surface. They are immediately bound to the outer layer of the cell wall or to only short filaments. On isolated cell walls the T protein was detected only on the outer surface.

CT Agglutination Tests
Antigens, Bacterial

Hines 09/518,165

*Bacterial Proteins

*Cell Wall

*Streptococcus pyogenes: UL, ultrastructure